

AIRSTAR[®]

DENTAL AIR SYSTEM
with Membrane Dryer



PRE-INSTALLATION GUIDE

Doctor:	_____
Address:	_____
Phone#:	_____
Dealer:	_____
Dealer Address:	_____


All installations must conform to local codes!

This AirStar Model is being installed:

(AS CHECKED)

- AS10 AS12 AS21 AS22
 AS30 AS40 AS50 AS70 AS100

Equipment Room Layout




Type		Style
AS10	5-20R	Green Dot
AS21	NEMA *	
AS12	6-20R NEMA *	Green Dot
AS22		
AS30		
AS40		
AS50	Hard Wired **	See Figure 3 for Connection Details.
AS70		
AS100		

* Use Hospital Grade Receptacle Only
 ** Only AS100 Models use 3-Phase Power.
 Disconnect Needed when Servicing.



Control Cable 18 Gauge
 4 Conductor Connect to
 Pressure Switch.
 See 24 Volt Circuit
 Connections.

Notes:

- Remote Air Intake Kit Manifold** - Refer to the table below listing the kit part number corresponding to the AirStar model.
- Remote Air Intake Drip Leg & Valve** - A drip leg with drain valve must be installed at lower end of the remote air intake pipe to collect condensation during operation. Attach a drain tube to the drip leg valve to allow drainage into floor drain/sink.
- Air Pipe** - Pipe can draw air from Inside or Outside. Drawing air from cooler and drier environments is beneficial to system performance. If drawing from Inside, the pipe end should still be protected from ingress of water and animals as appropriate. Be aware that noise can emit from the pipe during compressor operation; place end accordingly.
- Membrane Dryer Drain Valve** - Install the Drain Tube found in either accessory kit P/N 87133 or P/N 87134 used to drain moisture collected in the Membrane Dryer.



AIR PIPE
 2-Inch Pipe for Air Intake.
 Must be protected from rain and animals See note 3.



Remote Air Intake Kits:

Model	AS10	AS12	AS21	AS22	AS30	AS40	AS50	AS70	AS100
Part No.	85491	85491	85492	85492	85492	87361	85493	85494	85511

Air System Plumbing Connection

- 3/8" FNPT Shut-off valve and a 6 ft. pressure hose (supplied)
- Air distribution piping for all models - 1/2", type "L" or type "K" copper
- If pipe volume is too great, more than 235 in³ or more than 100 ft. of 1/2" diameter pipe, a pressure regulator should be installed between the main tank and the distribution piping and set to 80 psi.

Service Clearance

- Allow 12" on all sides for all units

Ambient Temperature

- Must not exceed 105°F
- Must be above 41° F




MEDICAL ELECTRICAL EQUIPMENT
 WITH RESPECT TO ELECTRICAL SHOCK, FIRE, MECHANICAL AND OTHER SPECIFIED
 HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES60601-1: A1:2012, C1:2009/(R)2012
 and A2:2010/(R)2012, CSA CAN/CSA-C22.2 NO.60601.14




Air Techniques, Inc.
 1295 Walt Whitman Road
 Melville, New York, USA 11747-3062

Electrical Connections




Notes:


1. Use 18 Gauge, 4 conductor, interconnect cable between the AirStar unit and Remote Switch Panel. When any 24V circuit exceeds 150 feet, use #16 AWG.
2. As shown, 3 conductors of the 24V circuit cable from each compressor connect via the user-supplied interconnect cable.
3. The fourth conductor of interconnect cable to be used for future equipment options/enhancements.
4. Leave with factory connection, without a 24V Switch, or connect the associated interconnect cable directly to remote switch. Please note that one switch is used for each compressor.



24V Circuit Connections



AS50 and AS70 Electrical Connections



AS100 3-Phase Electrical Connections

AirStar Electrical Connection

Site Requirements

	AS10	AS12	AS21	AS22	AS30	AS40	AS50	AS70	AS100
Note 1: Install a buck or boost transformer if actual facility service is above or below the supply voltage fluctuation of $\pm 10\%$ of nominal voltage ratings listed, install a buck/boost transformer with the corresponding part number as shown below.									
Nominal Supply Voltage (VAC, see note)	120	220	120	220	220	220	220	220	3~220
Maximum Current (Amps)	8	4	15	8	8	12	16	24	20
Minimum Panel Breaker Rating (Amps)	20	10	20	20	20	20	20	40	30
Minimum Wire Size (AWG)	12	16	12	12	12	12	12	8	10
Maximum Watts	864	792	1,620	1,584	1,584	2,376	3,168	4,752	3,960
BTU/hr @ 50%	1,473	1,350	2,201	2,701	2,701	4,051	5,401	8,201	6,752
BTU/hr @ 75%	2,210	2,026	4,143	4,051	4,051	6,077	8,102	12,153	N/A
Note 2: The 220 volt 3-phase, 4.6 KVa transformer is not available from Air Techniques.									
Buck/Boost Transformer Air Techniques Part No. and Size	67006 2.5 KVa	67002 3.9 KVa	67006 2.5 KVa	67002 3.9 KVa	67002 3.9 KVa	67002 3.9 KVa	67002 3.9 KVa	67000-1 7.8 KVa	4.6 KVa N/A

All AirStar compressors comply with NFPA 99C level 3 requirements

Product Specifications - Dimensions

Description \ Model	AS10	AS12	AS21	AS22	AS30	AS40	AS50	AS70	AS100	
Input Voltage Phase	1	1	1	1	1	1	1	1	3	
System Power HP (kW)	0.75 (0.56)	0.75 (0.56)	1.5 (1.1)	1.5 (1.1)	1.5 (1.1)	2.25 (1.68)	3.0 (2.2)	4.5 (3.4)	4.8 (3.6)	
Nominal Supply Voltage VAC	120	220	120	220	220	220	220	220	220	
Frequency Hz	60	60	60	60	60	60	60	60	60	
Maximum Current Amps	8	4	15	8	8	12	16	24	20	
Maximum Simultaneous Air Users	2	2	3	3	4	6	8	12	15	
System Output Flow Rate at 80 PSI CFM	2.5	2.5	5.0	5.0	5.0	7.5	10.0	15.0	20.0	
Maximum Pump-up Time 0-115 PSI (M:SS)	2:55	2:55	3:10	3:10	3:10	1:40	2:50	2:40	2:30	
Maximum Recovery Time 85-115 PSI (M:SS)	0:48	0:48	0:47	0:47	0:47	0:34	0:42	0:40	0:30	
Tank Size US Gal. (ft³)	6 (0.8)	6 (0.8)	12 (1.6)	12 (1.6)	12 (1.6)	12 (1.6)	20 (2.7)	30 (4.0)	30 (4.0)	
Shipping Weight lbs. (Approximate)										
No Sound Cover	170	170	200	200	240	255	290	430	450	
With Sound Cover	215	215	240	240	285	300	335	N/A	N/A	
Dimensions in. (See note.)										
No Sound Cover	H	28.50	28.50	30.50	30.50	30.50	30.50	33.50	35.00	34.50
	W	25.00	25.00	29.00	29.00	29.00	32.50	35.50	47.75	46.00
	D	19.75	19.75	20.00	20.00	20.00	20.00	20.50	21.75	24.50
With Sound Cover	H	30.00	30.00	32.00	32.00	32.00	32.25	33.50	36.00	N/A
	W	25.00	25.00	31.00	31.00	31.00	33.25	36.50	51.00	N/A
	D	22.50	22.50	22.25	22.25	22.25	22.50	22.75	29.50	N/A

Note: Height measured without leveling feet for all units with or without sound cover.

www.airtechniques.com



Corporate Headquarters
1295 Walt Whitman Road Melville, New York 11747- 3062
Phone: 800-247-8324 Fax: 888-247-8481

