PRO VECTA®

Panoramic & Cephalometric X-Ray Systems

Installation Planning Guide

Installation Prerequisites

Site Requirements

🔲 Closed, dry space n	ot made for anothe	r purpose (e.g. boiler	room or wet cell).
-----------------------	--------------------	------------------------	--------------------

- ☐ Have minimum footprint measuring 50 inches wide by 52 inches deep for ProVecta S-Pan and 80 inches wide by 52 inches deep for ProVecta S-Pan CEPH.
- ☐ The ceiling should be at least 8-foot high for both models.
- The flooring should be stable, level and capable of supporting 300 pounds equipment weight; 410 pounds with optional base plate assembly that adds 110 pounds.
- ☐ Have no large interference fields (e.g. strong magnetic fields) present, that can interfere with unit operation.
- ☐ Have dedicated 220V, 20 AMP, single phase 50/60 Hz circuit line prepared for connection.
- Environmentally controlled to temperatures between 50 and 95 °F with a relative humidity range of 30 to 75 %.

Recommended Computer System Requirements

Processor/CPU ≥ Intel Pentium IV compatible, Display Adapter 1

1.4 GHz

200 GB

Display Adapter 1280 x 1024 32bit color video display adapter (True

color) 128MB or greater

Direct3D®-capable workstation-class graphics card

RAM 4 GB or greater Operating System Microsoft Windows 8.1 (not Windows RT)

Microsoft Windows 10 (Pro or higher)
Microsoft Windows Server 2016
Microsoft Windows Server 2019

Network Interface Gigabit Ethernet adapter
Slots 1 PCI Express x 1
Optical Drive SuperMulti DVD Drive

Hard Disk Drive



Provecta S-Pan Model



Provecta S-Pan CEPH Model


www.airtechniques.com









ProVecta S-Pan Panoramic X-Ray Dimensions






Depth	Width	Height
В	Α	E
44.88 to 45.67	39 in.	27.56 in.
(114 to 116 cm)	(99 cm)	(70 cm)
IMPORTANT Measurements provide minimum clearances for operation of the Rotating Unit.	С	F
	21.25 in.	89.88 in.
	(54 cm)	(228 cm)
	D	G
	17.75 in.	62.48 in.
	(45 cm)	(159 cm)
		Н
		49.29 in.
		(125 cm)

Wall & Floor Preparation (Wall mount requirements are same for both S-Pan and CEPH units.)




Optional Base Plate


(Does not require wall mount)



ProVecta S-Pan Cephalometric X-Ray Dimensions



Height Values				
Α	D			
86.38 in.	89.88 in.			
(219.4 cm)	(228 cm)			
В	E			
58.82 in.	101.69 in.			
(149 cm)	(258 cm)			
С	F			
27.56 in.	62.32 in.			
(70 cm)	(158 cm)			



	Depth	Width		
	а	С	е	
۱,	44.88 to 45.67 in.	55.04 in.	21.25 in.	
	(114 to 116 cm)	(140 cm)	(54 cm)	
	b	d	f	
	31 in.	76.29 in.	39 in.	
	(79 cm)	(194 cm)	(99 cm)	



All measurements provide minimum clearance distances for operation of the Rotating Unit.