

# VistaPano S 2.0

Panoramic & Cephalometric X-Ray Systems

## Installation Planning Guide

### Installation Prerequisites

#### Site Requirements

- Closed, dry space not made for another purpose (e.g. boiler room or wet cell).
- Have minimum footprint measuring 50 inches wide by 52 inches deep for VistaPano S 2.0 and 80 inches wide by 52 inches deep for VistaPano S Ceph 2.0.
- The ceiling should be at least 8-foot high for both models.
- The flooring should be stable, level and capable of supporting 300 pounds equipment weight; 410 pounds with optional floor stand that adds 110 pounds.
- Have no large interference fields (e.g. strong magnetic fields) present, that can interfere with unit operation.
- Have dedicated 220 V, 20 AMP, single phase 50/60 Hz circuit line prepared for connection.
- Environmentally controlled to temperatures between 50 and 95 °F with a relative humidity range of 30 to 75 %.

#### Recommended Computer System Requirements

<b>Processor/CPU</b>	≥ Intel Core i5	<b>Display Adapter</b>	1280 x 1024 32-bit color video display adapter (True color) 128MB or greater
<b>RAM</b>	16 GB or greater	<b>64-bit Operating System</b>	Microsoft Windows 10 (Pro or higher) Microsoft Windows 11 (Pro or higher) Microsoft Windows Server 2019 Microsoft Windows Server 2022 Microsoft Windows Server 2025 (VistaSoft 4.0.11 or higher)
<b>Hard Disk Drive</b>	200 GB		
<b>Network Interface</b>	Gigabit Ethernet adapter		
<b>Slots</b>	1 PCI Express x 1		

VistaPano S 2.0



VistaPano S Ceph 2.0



[www.airtechniques.com](http://www.airtechniques.com)





**Distributed by:**  
Air Techniques, Inc.  
1295 Walt Whitman Road  
Melville, NY 11747-3062, U.S.A.  
Phone: +1-800-247-8324  
Fax: +1-888-247-8481  
[www.airtechniques.com](http://www.airtechniques.com)



**Manufactured by:**  
VATECH Co. Ltd.  
13, Samsung 1-ro 2-gil  
Hwaseong-si, Gyeonggi-do, 18449  
Korea  
Phone: +82 1588-9510  
[www.vatech.com](http://www.vatech.com)





## VistaPano S 2.0 Dimensions




Depth	Width	Height
<b>B</b> 45.08 to 45.67 (114.5 to 116 cm)	<b>A</b> 38.98 in. (99 cm)	<b>E</b> 27.56 in. (70 cm)
 <b>IMPORTANT</b> Measurements provide minimum clearances for operation of the Rotating Unit.	<b>C</b> 20.55 in. (52.2 cm)	<b>F</b> 90.04 in. (228.7 cm)
	<b>D</b> 18.43 in. (46.8 cm)	<b>G</b> 62.48 in. (158.7 cm)
		<b>H</b> 49.29 in. (125.2 cm)

## Wall & Floor Preparation

(Wall mount requirements are the same for both units.)




Use the provided aligning plate as a template for the 4 mounting holes in the floor.





## Optional Floor Stand

(Adds 1.1 in. (2.8 cm) to all heights)



## VistaPano S Ceph 2.0 Dimensions



Height	
<b>A</b> 89.21 in. (226.6 cm)	<b>D</b> 90.04 in. (228.7 cm)
<b>B</b> 61.65 in. (156.6 cm)	<b>E</b> 62.48 in. (158.7 cm)
<b>C</b> 27.56 in. (70 cm)	

Depth	Width	
<b>a</b> 45.08 to 45.67 in. (114.5 to 116 cm)	<b>c</b> 54.25 in. (137.8 cm)	<b>e</b> 20.55 in. (52.2 cm)
<b>b</b> 28.93 to 29.53 in. (73.5 to 75 cm)	<b>d</b> 74.80 in. (190 cm)	<b>f</b> 38.98 in. (99 cm)



Measurements provide minimum clearances for operation of the Rotating Unit.