## SepaStar amalgam separator

The following information supplements the installation and operating instructions of the unit.

As a general rule, the installation and operating instructions of the device must also be observed. These instructions include important information such as safety instructions and information on the setup, electrical connections, disinfection process, cleaning process etc.

The current version of the installation and operating instructions is available in the download Center:

 Ambient temperature

 Upper humidity limit

 Only qualified specialists or persons trained by Air Techniques may install, connect, and commission the unit.

### Scope of delivery

The following items are included in the scope of delivery (possible variations may apply due to country-specific requirements and/or import regulations):

**SepaStar amalgam separator . . . . . A1700**

- Amalgam separator
- Set of connection fittings
- Flow restrictor set
- Hose ø 20 mm
- Display panel
- Cable for display panel, 1 m
- Quick installation guide

### Installation

 Prior to working on the unit or in case of danger, disconnect it from the mains.

 Separation from mains power supply:  
The device does not have a main switch and must be disconnected from the power supply via an external isolating switch.

### Important information for installation


**Please read and apply the following information. This is an essential part of ensuring trouble-free operation.**

- At the water inlet of the SepaStar install a suitable flow restrictor that is appropriate for the installation conditions.  
For information about choosing the correct flow restrictor, please refer to the quick installation guide or the installation and operating instructions.
- Install the SepaStar as closely as possible to the Mojave Tank.
- The connecting hoses should be installed to be as short as possible. If necessary, shorten the hoses to the required length.
- Route the hoses with a continuous downward gradient towards the SepaStar. Avoid the formation of a siphon.
- The display must be mounted in such a way that the surgery personnel can see the signals at all times.


## Connecting the external start and display panel to the SepaStar

1  
Undo the four screws.

2  
Remove the cover towards the front.

### Cable entry and strain relief


1 Cable entry to the electronics  
2 Strain relief

### External start connection

1  
Guide the connection cable with the connector for external starting through the conduit to the electronics on the housing of the SepaStar.

2  
Connect the plug to connector X5.

3  
Fix the connection cable with strain relief.

### Connecting the display panel to the SepaStar

1  
Connect the cable with RJ-45 connector to the SepaStar electronics (X8) and plug it into the RJ-45 socket for direct connection to the display panel.

2  
Fix all cables with strain relief.

3  
Close the electronics cover of the SepaStar properly.

### Connecting the display panel

**i** The display panel is used to indicate messages acoustically and visually (via LEDs).

If the amalgam separator is set up at a larger distance from the treatment chairs (e.g. in the basement), the display panel must be installed appropriately such that the status messages of the amalgam separator can be monitored at all times.


**i** Installation, for example, on the wall next to or near low-voltage remote switches in the office.

**Important:** The LEDs on the display panel must be visible at all times.

**i** There must be a direct line connecting the RJ-45 socket on the unit and the RJ-45 socket on the display panel. Do not toggle network units (e. g. switch or router).

Pay attention to the resistance of the network cable between the RJ-45 sockets. The maximum length should not exceed 50 m.

1  
Connect display panel and RJ-45 socket using the ISDN cable supplied.

### Electrical connection

**WARNING**  
**Electric shock**

➤ The device may only be connected to a supply system with an earthed power outlet.

1  
Check conductivity setting (jumper position X10) and adjust as necessary.

2  
Establish the electrical connection to the supply network.


**i** Separation from mains power supply:  
The unit has no main power switch. An external mains disconnect must be installed upstream of the unit.

3  
Connect the display panel.

4  
Connect external start (optional, recommended for systems with separating tank).

5  
Connect external signal (optional).

6  
Connect the network (when monitoring via the network).

- X1 Mains AC power connection
- X2 Motor connection
- X3 Sensor system connection
- X4.1 Emergency start sensor
- X4.2 Reference sensor
- X4.3 Water start sensor
- X5 External starting (optional input, safety extra-low voltage 12-24 V, AC/DC)
- X6 External signal (switching capacity max. 24V, 120mA, AC/DC)
- X7 Relief valve connection
- X8 Display panel connection (RJ45 connector)
- X8a Display panel connection (6-pin connector)
- X10 Sensor sensitivity, conductivity value 100/200  $\mu\text{S}$
- X11 Network connection 100 Mbit
- X12 Diagnostic connector
- X14 Micro SD card holder
- F1 Fuse, brake, T 3.15 AH (IEC 60127-2)
- F2 Fuse, T 4.0 AH (IEC 60127-2)
- F3 Fuse, T 4.0 AH (IEC 60127-2)

**Network connection**

**i** All connected IT units must demonstrably comply with the current edition of IEC 62368.

**Purpose of the network connection**


The network connection is used to exchange information or control signals between the unit and a software installed on a computer, in order to, e. g.:

- Display parameters
- Select operating modes
- Indicate messages and error situations
- Change device settings
- Activate test functions
- Transmit data for archiving
- Provide documents concerning the devices

**Connecting the unit to the network**

**i** During initial installation, a router or server with DHCP is recommended so the unit is detected in the network.


- 1 Remove the cover from the electronics.
- 2 Plug the network cable into the electronics and into a network socket.
- 3 Attach the network cable to the device.
- 4 Create a connection to the network in the practice with the network cable.

**Installing the flow restrictor**

**i** To maintain a separation rate of at least 99%, the maximum flow rate of the liquid supplied must be reduced to 9 l/min.

The amount of fluid that flows into the amalgam separator depends on the height difference between the intake on the amalgam separator and the maximum water level of the connected devices. To limit the amount of fluid, an appropriate flow restrictor must be installed at the intake of the amalgam separator. Refer to the table below for help with selecting the required flow restrictor.

Height difference	Flow restrictor	
Min – 17.8"	#1	

Height difference	Flow restrictor	
17.9" – 34.6"	#2	
34.7" – 59.6"	#3	
> 59.7"	Not supported	


**i** Device combinations with corresponding flow restrictor can be found in the attachment.

**i** Flow restrictor #2 is already pre-installed, but it may need to be replaced by a suitable flow restrictor.

**1** Select the required flow restrictor to suit the installation conditions.

**2** Insert the flow restrictor into the fitting of the ventilation.

**3** Push the hose bushing onto the fitting and secure with a fixing ring.

### Connecting and routing hoses

**i** Install the hoses so they are as short as possible and with a gravity drain, without a siphon. This is the only way to avoid deposits in the hoses and ensure proper operation.


**1** Position the SepaStar as close as possible to the tank module.

**2** Cut the hoses to the required length.

**3** Screw hose sleeves onto the hose ends.

**4**


Connect the hoses to the hose connectors and secure them with hose clamps.

- 1 Hose adapter
- 2 Hose sleeve
- 3 Hose clamp
- 4 Outlet hose  $\varnothing$  20 mm
- 5 Inlet hose  $\varnothing$  20 mm

**5**

Attach the tank connection for the outlet hose to the Mojave tank.

- 1 Tank connection

**6**

Check the hose length and shorten if necessary.

**7**

Tighten hose sleeves onto the hose ends.

**8**

Attach the hose to the tank connection and secure with a hose clamp.

**9**


Run the outlet hose from the SepaStar with an incline to the outlet and connect it there.

### Connecting the external start to the Mojave

**1**

Loosen the two screws on the MMC-M cover.

- 2  
Open the housing door.


- 3  
Guide the connection cable for the external start to the MMC-M and through the strain relief for Aux/Alarm of the Mojave.

- 4  
Connect the plug (X2) labeled “MMC-M” on the connection card into the plug connector (J18) of the MMC-M.

- 5  
Plug the two flying leads of the connection cable into the plug connector (X1) labeled “START” at the top of the connection card and screw them together.

- 6  
Secure the strain relief. Ensure that the cable is connected without mechanical tension.

- 7  
Close the electronics cover of the Mojave properly.

## Commissioning and first start-up

- i** In various countries medical devices and electrical equipment are subject to regular checks at set intervals. The owner must be instructed accordingly.

- 1  
Turn on the unit power switch or the main surgery switch.

- 2  
Carry out the electrical safety check in accordance with national and local regulations.

- 3  
Check that coarse filters are installed in the units upstream of the amalgam separator.

- 4  
Check that the ventilation is installed and the correct flow restrictor has been used.

- 5  
Carry out a functional test.

- 6  
Check the unit and connections for leak tightness.

- 7  
Check that electronic covers are correctly mounted.

## Measuring the conductivity of tap water

- i** We recommend measuring the conductivity of the tap water at the installation site and adjusting the sensitivity of the conductivity sensor.

The sensitivity of the conductivity sensor for the water start can be adjusted using jumper X10 on the control panel. Jumper X10 is set by default on delivery to “100”.

Recommended set value:

- For conductivity > 320  $\mu\text{S}$ : set jumper to position “200”.
- For conductivity  $\leq$  320  $\mu\text{S}$ : set jumper to position “100”.

### When the jumper is initially set to “200”:

If the unit does not start after a period of time despite water supply, reduce the sensitivity of the conductivity sensor (jumper position “100”) and check the unit starts safely. This can be caused by seasonal fluctuations in the water supply or contamination deposits on the sensor.

## Monitoring the unit via the network

The following requirements must be met in order to monitor the device on the computer:

- Device connected to the network
- Current monitoring software installed on the computer

- i** As the monitoring system of the device, the software must deliver acoustic signals. Audio output on the computer must be activated.

## Combining devices safely

- Safety and essential performance features are independent of the network. The device is designed appropriately for operation independent of a network. However, some of the functions are not available in this case.
- Faulty manual configuration can lead to significant network problems. The expert knowledge of a network administrator is required for configuration.
- The data connection utilizes part of the bandwidth of the network. Interactions with other medical devices cannot be completely excluded. Apply the IEC 80001-1 standard for risk assessment.
- The device is not suitable to be connected directly to the public internet.

## Network configuration

Various options are available for network configuration:

- ✓ Automatic configuration via DHCP (recommended).
- ✓ Automatic configuration via Auto-IP for direct connection of unit and computer.
- ✓ Manual configuration.

- 1  
Configure the network settings of the unit using the software or, if applicable, the touch screen.

Check the firewall and release the ports, if applicable.

---

**Network protocols and ports**

Port	Purpose	Service
45123 UDP, 45124 UDP	Device recognition and configuration	
1900 UDP	Service indicator	SSDP / UPnP
502 TCP, 8080 <sup>1)</sup> TCP, 2005 TCP	Unit data	
514 <sup>1)</sup> UDP	Event protocol data	Syslog
22 TCP, 23 TCP	Diagnosis	SSH, Telnet
123 UDP	Time	NTP
2006	Diagnosis	


<sup>1)</sup> The port can vary depending on the configuration.

---

## Usage


### Ready for use

1


### Amalgam collecting container is 95% full


1

2


3


### Amalgam collecting container is 100% full


1

2


3


### Fill level measurement fault

1


### No amalgam collecting container inserted

1


2


### Motor fault

1


### Emergency start sensor in overflow position

1

## Cleaning

After each patient: suction 0.5 l of water

### Daily cleaning

Daily in the evening after finishing treatment; twice daily with high loads.

Always use:


- Monarch Evacuation System Cleaner
- Unit care system, e. g. Monarch Dispenser System

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