


How to store and maintain your equipment if you need to temporarily close your practice due to COVID 19.



Equipment

- Vacuum Lines
- Mojave Dry Vacuum
- Mojave LT3
- AirStar / AirStar NEO
- VacStar / VacStar NEO
- STS Dry Vacuum
- Remote Water Control Valve

Vacuum Lines

- a. Begin by running Clean Stream through the operatory furthest from your utility room (click link to watch instructional video <https://youtu.be/rFO-Kn08JYk>).
- b. Move to the next furthest and repeat
- c. Continue this process until all of your operatories have had their lines cleaned.
- d. Open an HVE in the operatory furthest from the equipment room (the one that was flushed in the first step) and allow the vacuum to run for 15 minutes. This will ensure that excess effluent is flushed from the lines.
- e. You may now begin shutting down your equipment.

Mojave Dry Vacuum

- a. Leave the unit connected to main power. Unit should be turned off at the low voltage remote control switch.
*Note the Mojave has a maintenance cycle that prevents the impeller from seizing up. The unit has to remain powered for this to happen.
- b. Water supply (for tank wash-out) should be turned off. This will remove the possibility of leaks.

Mojave LT3

- a. Turn off low voltage remote control switch. Unit must remain connected to power supply for maintenance cycle *see above.

VacStar / VacStar NEO

- a. Turn off water supply. Turn off main breaker.
- b. Replace or rinse out the solids collector. While not necessary this will help to prevent any odors building up in the equipment room during your shut down.

AirStar / AirStar NEO

- a. Turn off remote switch or disconnect from main power.

STS Dry Vac

- a. Can shut off remote switch or disconnect from power. Please shut off water supply for tank washout.

Remote Water Control Valve

- a. Turn off valve from main water supply. Turn off remote switch.

Digital

- AT2000 XR
- PeriPro
- ScanX
- ProVecta S-Pan
- ProVecta HD

AT2000 XR

- a. All Chemistry should be drained. Unit should have a quarterly maintenance performed, i.e. clean rollers and tanks as prescribed. Water supply to wash tank should be shut off.

PeriPro

- a. See above. Perform quarterly maintenance.

ScanX

- a. Simply turn the unit off. Employ a dustcover. "This could be as simple as putting a pillow case or garbage can liner over the unit"

ProVecta S-Pan

- a. Turn off and unplug from wall. (Turn off circuit breaker if hardwired)

ProVecta HD


- a. Turn off and unplug from wall. (Turn off circuit breaker if hardwired)

Merchandise

- Monarch Line Cleaner

Monarch Line Cleaner

- a. Flush lines as prescribed.
- b. Rinse, dry and re-attach bottle (water bottles should not be stored full. This can promote growth of bio material)
- c. Run syringe and hand piece tubings until air come out. This will ensure dry lines.
- d. Remove bottle and store safely.



Questions?

Air Techniques' Technical
Support Experts

1- 800-AIR TECH(247-8324)

