Infection Control Practices for Dental Radiology



Using ScanX Phosphor Storage Plates with a Barrier Envelope



- 1. Completely slide the PSP into the barrier envelope.
- 2. After the exposure, in the event of heavy soiling, e.g. from blood, dry clean the barrier envelope. Disinfect the barrier envelope and protective gloves with a disinfection wipe or a spray disinfectant on a soft, lint-free cloth. Place the barrier envelope with the PSP on the disinfection wipe, allow the barrier envelope to fully dry.
- 3. Slide the PSP out of its barrier envelope downwards into the device until the PSP is automatically drawn in.







Monarch Surface Disinfectant Wipes

- Kills TB, MRSA, HIV, HBV, Polio and 20 other pathogens in just 1 minute.
- Ethanol-based & bleach-free.
- Low toxicity rating.

IS A PROUD MEMBER OF







SCAN>X[®] Swift View



The digital radiography system with a touch screen, for all intraoral formats.

Benefits at a glance:

- High resolution image quality for easier diagnosis
- Intuitive, efficient and time-saving digitization of Phosphor Storage Plates (PSPs) for all intraoral formats (sizes 0 - 4)
- Handling is as easy as for analog film
- PSPs are simple to position and sit comfortably in the patient's mouth
- High resolution touchscreen
- PC connection via WIFI/LAN





Size 2





The ScanX Swift View enables the digitization of all Intraoral PSP formats. Handling is as easy as for analog film.





