

Air Techniques, Inc. Introduces Polaris™ Intraoral Camera

Date: 01/01/2009



MELVILLE, NEW YORK. Air Techniques, Inc. presents Polaris, a new generation of intraoral cameras that combines exceptional image quality and ergonomic excellence. The camera's patent-pending 330 degree swivel head enables you to capture all angles with ease and the rounded, ergonomic design makes Polaris easy to hold and operate. Polaris' optical glass lens has tremendous depth of field. All objects from 6 mm to 40 mm are automatically in sharp focus. The sharpness and true-to-life soft tissue color make Polaris a valuable diagnostic tool and a fantastic aid to patient presentation. Two models are available to suit your needs. Polaris USB provides power and data transfer in a single computer connection. And, Polaris Universal offers multiple connectivity options including USB 2.0, S-Video, and Composite Video Outputs.

For more information, visit the [Polaris product page](#) or call 1-800-AIR-TECH.

Founded in 1962, Air Techniques is a leading manufacturer of quality equipment for the dental, veterinary, medical, and security industries.
