

Name	Software Version	USB / EPPApproved	Windows Compatibility			Contact Info	DICOM Compatible		Specialty
			2000	XP	Vista**		Yes	No	
Apteryx, Inc.	Version 3.9 or later	X	Yes	Yes	Yes	Phone: (877) 278 3799 Email: info@apteryxware.com http://www.apteryxware.com		X	All
ADSTRA Imaging	ADSTRA Imaging Version 3.9	X	Yes	Yes	Yes	Phone: 1-800-450-0220 Email: info@adstra.com http://www.adstra.com		X	All
CaptureLinkK	CaptureLink Version 10+ .Net Version 2+	X	Yes	Yes	Yes	Phone: 1-800-668-5558 Email: custserv@henryschein.ca	X		All
Dentrix	Image Version 3.1.53 / 4.x	X	Yes	Yes	Yes	Phone:(800) 735-5518 Email: support@dentrix.com http://www.dentrix.com	X		All
Dolphin	Dolphin Versions 9.0/10.0	X	Yes	Yes	No	Phone: (800) 548-7241 Email: sales@dolphinimaging.com http://dolphinimaging.com		X	All
Dexis	Dexis Version 8.0.0/8.0.6	X	Yes	Yes	No	Phone: (888) 883-3947 http://www.dexis.com		X	All
IMAGIN Systems Corp.	Imagin SystemsDentalEye3 Version 3.0.471	X	Yes	Yes	No	Phone: (888) 478-4624Email: webmaster@imaginsystems.com http://www.imaginsystems.com	X		All
Patterson EagleSoft	EagleSoft Version 11.0/12.0/13.0/14.0	X	Yes	Yes	Yes	Phone: (800) 294-8504 Email: sales@eaglesoft.net http://www.eaglesoft.net	X		All
Image FX	Image FX Version 4.0 / 5.0	X	Yes	Yes	No	Phone: 1-800-572-1211 http://www.scican.com		X	All
Image XL	Version 3.9 or later	X	Yes	Yes	Yes	Phone: (800) 328-2925 Email: xldentinfo@xldent.com http://www.xldent.com	X		All
MediaDent	Mediadent Version 4.5 / 5.0	X	Yes	Yes	No	Phone: 877-770-8514 Email: sales@multimediodental.org http://www.mediadentusa.com	X		All
Medicor Imaging	MiPACS Dental Imaging Enterprise. Plug-in Ver. 1.2 or 2	X	Yes	Yes	No	Phone: 704-227-2624 Email: sales@medicorimaging.com http://www.mediacorimaging.com	X		Government/ Institutional
Mogo Inc.	Mogo 10.0 / 11.01	X	Yes	Yes	Yes	Phone: (800) 944-6646Email: mogo@mogo.com http://www.mogo.com		X	All
Patient Gallery	Version 1.45 or later	X	N/A	Yes	N/A	Phone: (866) 917-9345 http://www.patientgallery.com			
PracticeWorks* (a Practice Works, Inc. Comp)	PracticeWorks DICOM Imaging Version 2.1 & 4.9 / Ver. 6.7	X	Yes	Yes	Yes	Phone: 1-800-944 6365 Email: dental@practiceworks.com http://www.practiceworks.com		X	All
SoftDent*	SoftDent Version 9.7 Ver. 12.0.8.118	X	Yes	Yes	Yes	Phone: (800) 433 2409 http://www.softdent.com		X	All
TigerView	TigerView Version 5.x thru 7.x	X	Yes	Yes	Yes	Phone: (800) 385 9593 Email: sales@tigerview.com http://www.tigerview.com	X		All
Tracker	ChairSide Imaging 5	X	Yes	Yes	Yes	Phone: (800) 92-BRIDGE (922-7434)Email: sales@bridgenetwork.com http://www.bridge-network.com		X	All
Vipersoft (Integra Medical)	Vipersoft Version 3.1.53 or 4.0	X	Yes	Yes	Yes	Phone: (877) Viper GoEmail: info@vipersoft.com http://www.vipersoft.com (877) 847-3746		X	All
VISIX	Version 2x	X	Yes	Yes	Yes	Phone: 1 866 90-VISIX (84749) http://visixsoftware.com		X	All
XDR		X	Yes	Yes	Yes	Phone: (888) 937-9729 http://www.xdr.radiology.com			All

* Please consult your imaging software vendor for information about bridging capabilities

** Vista compatibility refers to the imaging software listed. ScanX has different compatibility restrictions. Contact software vendor for specific computer and detailed operating system requirements