









Authorized Imaging Software Vendors for SLC USB & Acclaim Intraoral Camera as of 10/18/2004

Name	Software Version No.	USB Approved	Windows Compatibility		Contact Info
			2000	XP	
 ADSTRA SYSTEMS <small>Operating Systems for High Tech Dentistry™</small>	ADSTRA Imaging Version 3.9	TBD	Yes	Yes	Phone: 1-800-450-0220 Email: info@adstra.com http://www.adstra.com
CaptureLink	CaptureLink Version 8.0+ available only in Canada	X	Yes	Yes	Phone: 416-573-8772 Email: albadent@cogeco.ca http://hsa.ca
 DENTRIX®	Image Version 4.0.65* or 4.0.83.2**	X	Yes	Yes	Phone: (800) 735-5518 Email: support@dentrix.com http://www.dentrix.com
 Patterson EagleSoft <small>PATTERSON DENTAL SUPPLY, INC.</small>	EagleSoft Version 9.0/10.0	X	Yes	Yes	Phone: (800) 294-8504 Email: sales@eaglesoft.net http://www.eaglesoft.net
 IMAGE FX 4.0	Image FX Version 5.0	X	Yes	Yes	Phone: 1-800-572-1211 http://www.scican.com
 practiceworks  SoftDent <small>(a PracticeWorks, Inc. company)</small>	PracticeWorks DICOM Imaging Version 2.1 & 4.9 SoftDent Version 9.7	TBD	Yes	Yes	Phone: 1-800-944 6365 Email: dental@praciceworks.com http://www.praciceworks.com Phone: (800) 433 2409 http://www.softdent.com
 sirona <small>The Dental Company</small>	Sidexis Version XG	X	Yes	Yes	Phone: 1-800-659-5977 http://www.sidexis.com
 TigerView <small>(Telvere Systems)</small>	TigerView Version 5.21/5.22	X	Yes	Yes	Phone: (800) 385 9593 Email: sales@tigerview.com http://www.tigerview.com
 Vipersoft <small>(Integra Medical)</small>	Vipersoft Version 4.0	TBD	Yes	Yes	Phone: (877) Viper Go Email: info@vipersoft.com http://www.vipersoft.com
VixWin	VixWin 2000 Version 1.8	X	Yes	Yes	Phone: (800)-800-2888 www.qendexxray.com
 XV3®	Version 3.5	X	Yes	Yes	Phone: (877) 278 3799 Email: info@apteryxware.com http://www.apteryxware.com

Computer Requirements/Recommendations:

CPU Speed:	2 GHz Pentium-4
Operating System:	Windows 2000 Service Pack 4 or Window XP Service Pack 1 with Microsoft knowledge base KB822603 update
RAM:	256 MB
Hard Drive	40 GB
CD ROM Drive:	48X
Monitor 1024 x 786 Resolution):	SVGA 17" (1024x768 resolution) Contrast ratio 450:1; .22 dot pitch capability CRT is recommended for optimum viewing
Video Display Adapter:	32 MB RAM
USB Port:	USB 2.0 High Speed Port
Keyboard, Mouse:	Standard
Printer:	Optional

*When Acclaim Intraoral Camera is used with Dentrix software version 4.0.65, please use the SLC driver disk in place of the Acclaim driver disk (the above software version only allows users to select the SLC camera and not the Acclaim camera).

** When Acclaim Intraoral Camera is used with Dentrix software version 4.0.83.2, please use the Acclaim driver disk.