

Spectra™ Caries Detection Aid



Key Features

- FDA Cleared
- TWAIN Compliant
- Hand-held
- Non-invasive imaging instrument
- Visually aids in the detection of caries
- Operates using fluorescence technology
- LEDs project high-energy light onto the tooth surface causing cariogenic bacteria to fluoresce red and healthy enamel green.

Read the Dental Product Shopper review of Spectra by clicking the image below!



[Spectra Caries Detector](#)

[Click here to see FAQs.](#)

[Click here to request a product demonstration.](#)

Product Specifications

Part Number	D4000
Dimensions (LxWxD)	8.6"x0.9"x1.2"
Weight	2.1 oz.

Connection Type	USB
Current	0.5 Amps
Package Includes	Handpiece, software, disposable camera covers (qty. 50), spacers (qty. 10), handpiece holder (1); USB cable

Computer Requirements

Computer Recommendations

CPU	Pentium 4, 2 GHz or higher
Operating System	Windows XP Professional with Service Pack 3 or later for an Intel 32-bit processor; Microsoft Windows XP Professional 64-bit Edition with Service Pack 2 or later for an Intel 64-bit extended (x64) processor; Microsoft Windows Vista Business, Enterprise, or Ultimate with Service Pack 1 or later for an Intel 32-bit or an Intel 64-bit extended (x64) processor; Microsoft Windows 7 Professional, Enterprise, or Ultimate for an Intel 32-bit or an Intel 64-bit extended (x64) processor.
RAM	256 MB
CD ROM Drive	24X speed
Monitor	SVGA 17", 800 x 600 or higher resolution or higher
Video Display Adapter	16 MB Video card supporting 800 x 600 pixels with 32-bit color
Keyboard/Mouse	Standard keyboard and mouse
USB Port	Must be USB 2.0 high speed, high power port
Printer	Optional

Video Links

Optional Accessories

- Disposable Camera Sheaths
- Spacers

Warranty Information

3 Year Warranty

This Air Techniques dental product is warranted to be free from defects in material and workmanship for a period of 3 years from the date of installation by authorized Air Techniques' dealer service personnel.

Downloads



[Spectra Brochure](#)



[Spectra Brochure \(Spanish\)](#)



[Spectra Caries Detection Aid - Driver \(.zip 90MB\)](#)



[Spectra Caries Detection Aid - Operator's Manual](#)



[Spectra Driver Instructions](#)



[Spectra Quick Start Guide](#)



[Spectra Sell Sheet](#)