

Choose Acadia™ Amalgam Separator, the All-Around Safe and Effective choice for Amalgam Collection.

Air Techniques is the industry leader for nearly five decades in powering dental practices with first-rate air and evacuation systems. At Air Techniques, we understand how to create an amalgam separator that works in-line with your evacuation system. Acadia is ISO 11143 certified and exceeds standards, ensuring compliance with all US and Canadian regulations.

Easy, no-hassle installation, regardless of space or plumbing configuration.

Acadia's small footprint allows for a snug fit and can be installed free-standing or wall-mounted. No additional electricity, wires or plumbing is needed and rotatable connections ensure your plumbing stays intact.

Smart design for your practice.

Available in a single, smartly designed model, Acadia supports up to ten operators. A translucent container enables you to view Acadia's cyclonic separation of fluids – including a clear look at the filter. An inlet filter ensures a smooth and continuous filtration process; all particles larger than 5mm are filtered out, never clogging or causing loss of suction.

Simple and safe recycling procedures.

Recommended to be changed annually, the four stage filter easily twists on and off making replacement simple (no tedious clips to fit into place). The visual fill line lets you know when your filter should be replaced. With each Recycle Kit sold, Air Techniques provides a prepaid shipping label and the necessary packaging for the filter to be sent to a certified recycler, along with a recycling certificate issued to your practice.



Acadia Amalgam Separator
PN A1250

**AIR
TECHNIQUES**
INC.

www.airtechniques.com

Acadia™ Alert keeps you out of the utility room.



Acadia Alert
PN A1330

Acadia Alert automatically monitors your system and activates an audible alarm along with visual LED indicators when the filter reaches approximately 90% capacity. Acadia Alert can be installed in any operatory to give you real-time notification. This system requires no additional electricity and is powered by your vacuum's 24V control circuit. Acadia Alert can also be integrated with the Vacuum Remote Switch and Vacuum Equalizer from Air Techniques.



Take the next step toward managing an eco-efficient practice with Acadia Amalgam Separator. To locate an authorized dealer near you, visit www.airtechniques.com.

Filter Change in 3 Easy Steps

Step 1: Twist filter off



Step 2: Place used filter in proper shipping materials



Step 3: Twist new filter on



Replacement Filter
PN A1290

PN	Description
A1250	Acadia Amalgam Separator
A1290	Acadia Filter Replacement Kit w/o recycling
A1330	Acadia Alert
A1390	Acadia Filter Replacement Kit w/ recycling
57850	CleanStream Cleaner

Technical Specifications	
Services	Up to 10
Liquid Flow Rate	1 liter per minute
Efficiency Rate	ISO 11143
Connection Type	1½" Hose barb adaptor with retaining clip
Dimensions (WxDxH)	21" x 15" x 27"
Installation	Free standing or wall-mounted

Keep your utility room at peak performance with CleanStream™.

CleanStream maintains your wet or dry vacuum system by preventing debris from settling in the lines. To keep your vacuum and overall system operating at peak performance, we recommend regular cleaning of your evacuation lines using CleanStream Cleaner. For prolonged life and uninterrupted workflow, choose CleanStream pH-balanced evacuation system cleaner.



CleanStream Cleaner
PN 57850

1295 Walt Whitman Road, Melville, New York 11747-3062
516-433-7676 • Fax: 516-433-7684



www.airtechniques.com