



CLEANSTREAM™ EVACUATION SYSTEM CLEANER

Your dental vacuum is an integral part of keeping your practice running smoothly and with CleanStream Evacuation System Cleaner you can ensure you are getting the most out of your investment. Protect your equipment, patients and the environment with CleanStream.

Why Use an Evacuation System Cleaner?

Organic and inorganic debris promotes the growth of disease and can cause system failure, which is why every vacuum should be cleaned daily to maintain peak performance. Removing these excess particles protects as well as extends the life of your vacuum pump.

Unlike other Evacuation System Cleaners that use citric acid or chlorine based cleaners, CleanStream keeps your vacuum system running smoothly and provides you with peace of mind that your vacuum is not only efficient but also safe for the environment.

Now in a New and Improved Bottle

- Compact and ergonomic shape for easier handling and dosing
- Maximum safety with new tamper proof closure
- Integrated measuring cap



Works with any Wet or Dry Vacuum and Amalgam Separators.



AIR TECHNIQUES INC.



CLEANSTREAM™ EVACUATION SYSTEM CLEANER



Good For You, Good For the Environment

- Special non-foaming agents reliably remove and prevent debris from settling in evacuation lines.
- Specially designed dispenser allows the cleaning of two operatories with one bottle - no measuring or monitoring needed!
- Will not cause leaching of mercury from amalgam.
- Environmentally friendly and completely biodegradable.

Part Number	Description
57850	CleanStream Vacuum System Cleaner 2.5 Liter bottle (125 applications)
57665	CleanStream Dispenser System Includes Dispenser System with fittings and instructions
57660	CleanStream Starter Kit Includes 1 dispenser system with all accessories and 6 applications of CleanStream Evacuation System Cleaner



Easy to Use:

Fill CleanStream Dispenser with Water



Add CleanStream concentrate



Connect hoses to automatically aspirate solution

Air Techniques is a leading manufacturer of dental equipment for nearly 50 years. We have earned our reputation through our commitment to ensuring customer satisfaction, and we are proud to provide superior quality to dealers and dental practitioners worldwide.

Air Techniques products are distributed through authorized dealers only.

Headquarters

1295 Walt Whitman Road, Melville, NY 11747-3070
Phone: 516-433-7676 □ **Fax:** 516-433-7684
Toll Free in the U.S. and Canada: 1-800-AIRTECH

Western Facility

291 Bonnie Lane, Suite 101, Corona, CA 92880-2804 □ **Phone:** 951-898-8555
Fax: 951-898-7646 □ **Toll Free in the U.S. and Canada:** 1-800-AIRTECH

AIR TECHNIQUES INC.

www.airtechniques.com

