

STS™ and AirStar® Enhancements



STS 6 w/ 8 gallon CAS



AirStar 30



New Membrane Dryer System
(Integration begins Spring 2009)

MELVILLE, NEW YORK. Air Techniques Inc. is proud to inform you of some exciting new enhancements to our existing Compressor and Vacuum product line. Air Techniques' air compressor and dry vacuum Systems are now manufactured with important engineering innovations designed to boost the overall quality and efficiency of these popular Air Techniques products

STS

- *Horizontal motor design allows for space saving flexibility
- *Nickel-plated housing and impeller reduces oxidation increasing efficiency and motor life
- *Larger bearing size increases reliability
- *CAS Washout Timer Control prevents debris build-up in the tank

AirStar

Quieter Operation

- * Improved intake filters
- * Redesigned Purge Muffler
- * NEW Rubber Shock Mounts
- * Flat frame reduces vibration

New Membrane Dryer System

- *Maximum dryness with triple-filtered air
- * Uninterrupted availability
- * Low pressure dew point
- *Compact shape for space-saving installation
- * AS 30 and 50 models available Spring 2009, complete line Summer 2009