

Monarch Cleanstream Evacuation System Cleaner



Air Techniques is North America's leading manufacturer of dental vacuum systems. Monarch CleanStream has been carefully formulated from years of experience and research to help you optimize system performance and keep your vacuum system running smoothly.

- Unique formulation removes and prevents debris build-up in evacuation lines maintaining peak performance and suction.
- Cleans and deodorizes evacuation lines for all wet and dry vacuum systems.
- Non-foaming, highly concentrated formula is very economical to use.
- Amalgam separator safe with a pH of 7-8.
- Non-corrosive, for daily use.

Monarch CleanStream is compliant with the EPA Rule BMP for Dental Amalgam Waste

The EPA passed its final rule specific to Best Management Practices for Dental Amalgam Waste on July 14, 2017. The final rule prohibits the use of bleach or chlorine-containing cleaners that may lead to the dissolution of solid mercury when cleaning chair-side traps and vacuum lines. "...vacuum lines that discharge amalgam process wastewater to a POTW [publicly owned treatment works] must not be cleaned with oxidizing or acidic cleaners, including but not limited to bleach, chlorine, iodine and peroxide that have a pH lower than 6 or greater than 8." (40 CFR 441.30(b)(2)) Monarch CleanStream complies with this rule with a pH of 7-8 (diluted) and does not contain any oxidizing agents.

This statement along with the complete rule, can be found on the EPA website:

<https://www.epa.gov/eg/dental-effluent-guidelines>



 **Monarch**

by

AIR
TECHNIQUES equipped for life®

For more information visit www.airtechniques.com or call our Monarch Support Team at **1-800-AIRTECH**

PN 9922-34 RevB