### Film Processor Maintenance

**IMAGING** | UTILITY ROOM | MERCHANDISE

YOU DEPEND ON RADIOGRAPHIC FILMS FOR DIAGNOSIS AND OFTEN EMPLOY THEM IN THE PRESENTATION OF TREATMENT PLANS OR AS A PATIENT EDUCATION AID.



# YOU REALIZE THE IMPORTANCE OF CONSISTENT, ARTIFACT FREE FILM PROCESSING.

- Regular maintenance of your processor with the right maintenance products ensures best quality radiographs and longest service life.
- Air Techniques' maintenance products are specially formulated for maximum cleaning efficiency and ease-of-use.
- Build-up residues are easily removed and staining is reduced, protecting films from roller deposits.
- All Air Techniques' maintenance products are environmentally safe; they do not contain any heavy metals.

#### FORMULA 2000 PLUS

The non-toxic and non-hazardous cleaning concentrate effectively removes deposits and stains from roller transports, processors and dip tanks.





## Film Processor Maintenance

# Maintenance for Roller Transport Processors with Formula 2000 PLUS and Spray 2000 PLUS

#### Quarterly Cleaning with Formula 2000 PLUS

- Reduces staining and prolongs the life of your processor.
- Non-toxic and non-hazardous.
- Effectively dissolves deposits, dirt, algae and residue.
- Deep cleaning of the Roller Transport Processors takes only 15 minutes.

#### Weekly Cleaning with Spray 2000 PLUS

- Removes residue and protect films from roller deposits.
- Available in a convenient, ready to use spray bottle.
- No mixing necessary just spray, sponge and rinse under running water.
- For use in A/T 2000 XR.



(Left): Formula 2000 Plus Cleanser Twin Pack, PN 43945 (Right): Spray 2000 Plus, PN 43965

# Maintenance for Peri-Pro Transports with Peri-Pro Transport Cleaning Kit

#### Peri-Pro Transport Cleaning Kit with Formula 2000 PLUS

- Reduces staining, decreases normal maintenance time and prolongs the life of your Peri-Pro.
- Easy and convenient way to care for your Peri-Pro.
- Specially designed to hold the transport with just the right amount of Formula 2000 PLUS and water.
- Tub is completely self-contained, portable and requires little counter space.
- Deep cleaning only takes 30 minutes.
- Can be used for all models of Peri-Pro.



Peri-Pro Transport Cleaning Kit, PN 43975





### www.airtechniques.com

#### **Corporate Headquarters**

1295 Walt Whitman Road | Melville, NY 11747-3062 Phone 800.247.8324 | Fax 888.247.8481

#### **Western Facility**

291 Bonnie Lane, Suite 101 | Corona, CA 92880-2804 Phone 800.247.8324 | Fax 951.898.7646

