



4+ ENZYMATIC CLEANER / SAFETY TRAY

For use in Ultrasonic Cleaners and as a Presoak

Superior cleaning and protection of your
valuable dental instruments.

 **Monarch.**

by


AIR
TECHNIQUES equipped for life®

Uniquely formulated with 5 active enzymes

Before



After



When challenged with a worst-case scenario, Monarch 4+ Enzymatic Ultrasonic Cleaner was able to successfully remove organic debris from all instruments.

Source:
*Dental Advisor Microbiology
Research Report,
Number 89 – October, 2019*

Instrument Protection

- It won't corrode or stain metals, even when carbide, stainless or aluminum instruments are mixed, resulting in protection of valuable instruments for longer life.

Surfactant Synergy

- Contains a blend of anionic and nonionic surfactants. These surfactants ensure that the instruments are clean and free of fragments that can redeposit during the cleaning process.

Efficient Process

- Liquid mixes instantly and dissolves completely without leaving any residue in the ultrasonic cleaner compared to tablets.
- Non-foaming with a fresh pleasant scent.

“Monarch 4+ enzymatic cleaner is the best instrument cleaner I've come across. It's easy to use, cleans the instruments the first time through, smells fresh, and isn't thick or gunky like some other products I've tried. It only takes ½ a cap full, so a little bit goes a long way.”

Karen Webster RDH, PHDH
Tri City Health Partnership
Dental Program Manager

Order Information

Order No.	Product Description
H6201	Monarch 4+ Enzymatic Cleaner, Bottle, individual 84.5 oz. bottle (makes 125 gallons); (4 bottles per shipper case)
H6310-N	Monarch 4+ Enzymatic Cleaner, Box, individual 50 x 1/3 oz. unit dose pouches, each makes 1-gal.(12 boxes of 50 pouches each per shipper case)
H6150	Monarch Safety Tray, 3 Quart Bath



CORPORATE HEADQUARTERS
1295 Walt Whitman Road
Melville, New York 11747- 3062
800-247-8324



Monarch is a registered trademark
of Air Techniques, Inc.
© Copyright 2021
PN 9922-606 RevA

