


P/N 53020
P/N 53170


P/N 53021
P/N 53173





V3M/V5M/V7M PUMP

FOR USE WITH
AIR TECHNIQUES
WATER SYSTEMS ONLY


WATER SWITCH



Use upper right ground stud with all solenoids and drain pump




OPTIONAL



To Single Color Switch
Connect by wire color


Use strain relief nearest to
respective terminal block.
DO NOT CROSS FEED WIRES

MASTER CONTROL PANEL SWITCH




P/N 53201-1 24V GREEN

See users manual for detailed wiring schematic



P/N 53202-1 6V GREEN/YELLOW

FLOAT SWITCH



FLOAT SWITCH



MT10-M

MOJAVE[®] MONITOR

MMC-M WIRING FOR MT10-M

P/N H5268M#B