


This document provides the connections necessary to install the MOJAVE DryStar Monitor System into a utility room that is properly prepared for the system to be installed.



Typical MOJAVE DryStar Monitor Single Pump Connection Overview

## **System Setup Connections**

Important:

Make sure to install the system in accordance with all local electrical and plumbing codes.

Make sure to efficiently use space by making connections as short and direct as possible to meet your particular site requirements.


Make sure that all hose connections are straight and secure without any sharp bends or kinks.


Since the vacuum hose is rigid, make sure not to stress connections especially at the pump inlet.


## MOJAVE Documentation.


Refer to the following Mojave documentation for more detailed information on the plumbing, electrical connections and site requirements as necessary.

- User's and Installation Manual, P/N H5545,
- Pre-Installation Guide, P/N H5546,
- CT20 Dual Tank Product Overview









## www.airtechniques.com

Corporate Headquarters

1295 Walt Whitman Road | Melville, New York 11747- 3062, USA Phone: 1-800-247-8324 | Fax: 1-888-247-8481







