



CLEANSTREAM & DISPENSER

Evacuation System Cleaner & Dispenser with **VORTEX** TECHNOLOGY

Optimizes suction performance keeping your vacuum running at its best.



Monarch.

by

AIR
TECHNIQUES equipped for life®

CleanStream is the ONLY evacuation line cleaner Air Techniques recommends using daily to ensure highest performance and longevity of your equipment.



Superior cleaning effectiveness

- A proprietary blend of surfactants and chelating agents, removes debris and prevents future build-up maintaining suction performance to ensure an uninterrupted workflow.
- Ideal cleaning results are achieved when CleanStream is used together with the Dispenser with Vortex technology.
- Works with both wet and dry vacuums.

Safe for you and your equipment

- Cleaning suction lines daily with Vortex action ensures optimal aerosol mitigation during procedures.
- Non-foaming solution protects vacuum for longer life.
- Non-corrosive & non-hazardous.
- Amalgam separator safe with a neutral pH.

Cost-Effective

- Only 20 ml of the concentrate is needed to clean 1 operator and deodorize lines.
- 125 cleanings per 84.5 oz. bottle.

“Since using CleanStream, our evacuation lines have been completely open, and we haven't had any clogged lines.”

Chris Bethel, RDH, PHDH

The CleanStream Dispenser with Vortex Technology Ensures Clean Vacuum Lines Throughout The Office



To see the full video, go to: <https://youtu.be/Z84Z5BZLgio>

Order No.	Product	Description
57660	Monarch CleanStream Starter Kit	Includes: 1 Dispenser and 1 bottle (33.8 oz.) of CleanStream
57850	Monarch CleanStream Evacuation System Cleaner	84.5 oz. Bottle with integrated measuring cap - 125 cleanings
57840	Monarch CleanStream Evacuation System Cleaner	33.8 oz. Bottle with integrated measuring cap – 50 cleanings
57083	Monarch CleanStream Evacuation System Cleaner	338 oz. Container – 500 cleanings
57665	CleanStream Dispenser with Vortex Technology	1 Dispenser with adapters for all suction tubings
H6210	Dispensing Pump (20 ml)	For CleanStream 338 oz. Container
H6209	Dispensing Aid (Canister Tap)	For CleanStream 338 oz. Container



CORPORATE HEADQUARTERS
1295 Walt Whitman Road
Melville, New York 11747- 3062
800-247-8324



Monarch is a registered trademark
of Air Techniques, Inc.
© Copyright 2023
PN 9922-619 Rev C

