SAFETY DATA SHEET 1. Product and Company Identification

	1. Product and Company Ident	incation
Product identifier	Monarch Disinfectant Wipes	
Other means of identification		
Synonyms	EPA Reg. No. 88494-2-89494	
Recommended use	Disinfectant	
Recommended restrictions	None known.	
Manufacturer information	Air Techniques, Inc. 1295 Walt Whitman Road Melville, NY 11747-3062 Phone: 516-214-5631 FAX: 516-433-2768	
Supplier	See above.	
	2. Hazards Identification	n
Physical hazards	Flammable liquids	Category 2
Health hazards	Acute toxicity, inhalation	Category 4
	Serious eye damage/eye irritation	Category 2A
	Specific target organ toxicity, single exposure	
	Specific target organ toxicity, single exposure	
Environmental hazards	Not classified.	
WHMIS 2015 defined hazards	Not classified	
Label elements		
Signal word	Danger	
Hazard statement	Highly flammable liquid and vapor. Causes ser respiratory irritation. May cause drowsiness or	rious eye irritation. Harmful if inhaled. May cause rdizziness.
Precautionary statement		
Prevention	Keep container tightly closed. Ground and bor	pen flames and other ignition sources. No smoking. nd container and receiving equipment. Use quipment. Use non-sparking tools. Take action to
Response	In case of fire: Use appropriate media to extinguish. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell.	
Storage	Store in a well-ventilated place. Keep cool. Ke	ep container tightly closed. Store locked up.
Disposal	Dispose of contents/container in accordance with local/regional/national/international regulations.	
WHMIS 2015: Health Hazard(s) not otherwise classified (HHNOC)	None known	
WHMIS 2015: Physical Hazard(s) not otherwise classified (PHNOC)	None known	
Hazard(s) not otherwise classified (HNOC)	None known.	

media

Although the product as a whole is in solid format, the product does not meet the OSHA HCS definition of a flammable solid as per Appendix B to 1910.1200 - Physical Hazard Criteria, section B.7.1 and B. 7.2 or Canada's Hazardous Products Regulations Subpart 7.

US: This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines.

This chemical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS), and for workplace labels of non-disinfectant chemicals. The hazard information required on the disinfectant label is reproduced below. The disinfectant label also includes other important information, including directions for use. PRECAUTIONARY STATEMENTS:

Hazards to humans and Domestic animals.

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating and drinking, chewing gum, using tobacco or using toilet.

Flammable. Keep away from heat or open flame.

3. Composition/Information on Ingredients

Mixture			
Chemical name	Common name and synonyms	CAS number	%
Ethanol		64-17-5	72.5
Didecyl dimethyl ammonium chloride		7173-51-5	0.33

All concentrations are in percent by weight unless ingredient is a gas. Gas concentrations are in percent by volume.

4. First Aid Measures		
Inhalation	IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell.	
Skin contact	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.	
Eye contact	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.	
Ingestion	Rinse mouth. Do not induce vomiting. Get medical attention if symptoms occur. Never give anything by mouth if victim is unconscious or is convulsing. Obtain medical attention.	
Most important symptoms/effects, acute and delayed	Symptoms may include stinging, tearing, redness, swelling, and blurred vision. Vapors have a narcotic effect and may cause headache, fatigue, dizziness and nausea.	
Indication of immediate medical attention and special treatment needed	Provide general supportive measures and treat symptomatically. In case of shortness of breath, give oxygen. Symptoms may be delayed.	
	Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center [NPIC] at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time [NPIC Web site: www.npic.orst.edu].	
General information	Ensure that medical personnel are aware of the material(s) involved and take precautions to protect themselves. Take off all contaminated clothing immediately. If you feel unwell, seek medical advice (show the label where possible). Show this safety data sheet to the doctor in attendance. Wash contaminated clothing before reuse. Keep away from sources of ignition. No smoking. Avoid contact with eyes and skin. Wear safety glasses with side shields. Protective gloves recommended based on infection control best practices. Keep out of reach of children.	
	5. Fire Fighting Measures	
Suitable extinguishing media	Dry chemical, CO2, water spray or regular foam.	
Unsuitable extinguishing	Do not use water jet as an extinguisher, as this will spread the fire.	

Specific hazards arising from the chemical	Vapors may form explosive mixtures with air. Vapors may travel considerable distance to a source of ignition and flash back. This product is a poor conductor of electricity and can become electrostatically charged. If sufficient charge is accumulated, ignition of flammable mixtures can occur. To reduce potential for static discharge, use proper bonding and grounding procedures. This liquid may accumulate static electricity when filling properly grounded containers. Static electricity accumulation may be significantly increased by the presence of small quantities of water or other contaminants. Material will float and may ignite on surface of water. During fire, gases hazardous to health may be formed. Firefighters should wear a self-contained breathing apparatus.	
Special protective equipment and precautions for firefighters	Firefighters should wear full protective clothing including self-contained breathing apparatus.	
Fire-fighting equipment/instructions	In case of fire and/or explosion do not breathe fumes. Move containers from fire area if you can do so without risk.	
Specific methods	Use standard firefighting procedures and consider the hazards of other involved materials.	
General fire hazards	Solid containing flammable liquid	
Hazardous combustion products	May include and are not limited to: Oxides of carbon.	
	6. Accidental Release Measures	
Personal precautions, protective equipment and emergency procedures	Immediately evacuate personnel to safe areas. Keep unnecessary personnel away. Keep out of low areas. Keep people away from and upwind of spill/leak. Eliminate all ignition sources (no smoking, flares, sparks, or flames in immediate area). Remove all possible sources of ignition in the surrounding area. Wear appropriate protective equipment and clothing during clean-up. Fully encapsulating, vapor protective clothing should be worn for spills and leaks with no fire. Do not touch damaged containers or spilled material unless wearing appropriate protective clothing. Avoid inhalation of vapors or mists. Avoid inhalation of vapors. Ventilate closed spaces before entering them. Use appropriate containment to avoid environmental contamination. Transfer by mechanical means such as vacuum truck to a salvage tank or other suitable container for recovery or safe	

Methods and materials for
containment and cleaning uppersonal protection, see section 8 of the SDS.Pick up and discard.
Eliminate all ignition sources (no smoking, flares, sparks, or flames in immediate area). Take
precautionary measures against static discharge. Use only non-sparking tools. Keep combustibles
(wood, paper, oil, etc.) away from spilled material.
Large Spills: Stop leak if you can do so without risk. Dike the spilled material, where this is
possible. Cover with plastic sheat to prevent spreading. Use a non-combustible material like

possible. Cover with plastic sheet to prevent spreading. Use a non-combustible material like vermiculite, sand or earth to soak up the product and place into a container for later disposal. Absorb in vermiculite, dry sand or earth and place into containers. Prevent entry into waterways, sewer, basements or confined areas. Following product recovery, flush area with water.

disposal. Local authorities should be advised if significant spillages cannot be contained. For

Small Spills: Wipe up with absorbent material (e.g. cloth, fleece). Clean surface thoroughly to remove residual contamination.

Never return spills to original containers for re-use.

Environmental precautions Avoid discharge into drains, water courses or onto the ground. Use appropriate containment to avoid environmental contamination.

7 Handling and Storage

	7. Handling and Storage
Precautions for safe handling	Vapors may form explosive mixtures with air. Do not handle, store or open near an open flame, sources of heat or sources of ignition. Protect material from direct sunlight. Do not smoke. Minimize fire risks from flammable and combustible materials (including combustible dust and static accumulating liquids) or dangerous reactions with incompatible materials. Handling operations that can promote accumulation of static charges include but are not limited to: mixing, filtering, pumping at high flow rates, splash filling, creating mists or sprays, tank and container filling, tank cleaning, sampling, gauging, switch loading, vacuum truck operations. Take precautionary measures against static discharges. All equipment used when handling the product must be grounded. Use non-sparking tools and explosion-proof equipment. Use only with adequate ventilation. Avoid breathing vapor. Avoid contact with eyes, skin and clothing. Avoid prolonged exposure. Use only outdoors or in a well-ventilated area. Wear appropriate personal protective equipment. Observe good industrial hygiene practices. Wash thoroughly after handling. Keep container tightly closed. When using do not eat or drink.
Conditions for safe storage, including any incompatibilities	Store locked up. Keep away from heat, sparks and open flame. Prevent electrostatic charge build-up by using common bonding and grounding techniques. Avoid spark promoters. Eliminate sources of ignition. Ground/bond container and equipment. These alone may be insufficient to remove static electricity. Store in a closed container away from incompatible materials. Store in a cool, dry place out of direct sunlight. Store in a well-ventilated place. Store away from incompatible materials (see Section 10 of the SDS). Keep out of reach of children. PESTICIDE STORAGE: Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

	8. Exposure Controls/Pe	ISONAL PROTECTION	
Occupational exposure limits			
•	ccupational Health & Safety Code, Scl		
Components	Туре	Value	
Ethanol (CAS 64-17-5)	TWA	1880 mg/m3	
		1000 ppm	
Canada. British Columbia Safety Regulation 296/97		s for Chemical Substances, Occupational Health and	
Components	Туре	Value	
Ethanol (CAS 64-17-5)	STEL	1000 ppm	
Canada. Manitoba OELs	Reg. 217/2006, The Workplace Safety	And Health Act)	
Components	Туре	Value	
Ethanol (CAS 64-17-5)	STEL	1000 ppm	
Canada, Ontario OFI s. ((Control of Exposure to Biological or C	hemical Agents)	
Components		Value	
Ethanol (CAS 64-17-5)	STEL	1000 ppm	
, , , , , , , , , , , , , , , , , , ,	Ainistry of Labor - Pogulation Posses	ting the Quality of the Work Environment)	
Canada. Quebec OELS. (I Components	Type	Value	
Ethanol (CAS 64-17-5)	TWA	1880 mg/m3	
		1000 ppm	
Components	ts for Air Contaminants (29 CFR 1910. Type	1000) Value	
Ethanol (CAS 64-17-5)	PEL	1900 mg/m3 1000 ppm	
US. ACGIH Threshold Lin		M. L.	
Components	Туре	Value	
Ethanol (CAS 64-17-5)	STEL	1000 ppm	
US. NIOSH: Pocket Guide	to Chemical Hazards		
Components	Туре	Value	
Ethanol (CAS 64-17-5)	TWA	1900 mg/m3	
		1000 ppm	
iological limit values	No biological exposure limits noted	for the ingredient(s).	
ppropriate engineering		0 air changes per hour) should be used. Ventilation rates	
ontrols		applicable, use process enclosures, local exhaust ventilation	
		ntain airborne levels below recommended exposure limits. If plished, maintain airborne levels to an acceptable level.	
dividual protection measure	es, such as personal protective equipr	· · · · · · · · · · · · · · · · · · ·	
Eye/face protection	Wear safety glasses with side shield		
Skin/Hand protection			
	on infection control best practices.	on test data. Protective gloves recommended based	
ther Respiratory	Where exposure guideline levels may be exceeded, use an approved NIOSH respirator.		
otection		· · · · ·	
ermal hazards	Not applicable.		
eneral hygiene	When using do not smoke. Always o	observe good personal hygiene measures, such as washing	
onsiderations		re eating, drinking, and/or smoking. Routinely wash work	
	clothing and protective equipment to	o remove contaminants.	
	9. Physical and Chem	ical Properties	
ppearance	-	•	
hysical state	Clear Liquid saturated on wipe		
orm	Liquid.		
Color	Liquid saturated on wipe		

Colorless

Not available.

Alcohol

Color

Odor

Odor threshold

рН	3.7 @ 77°F (liquid)		
Melting point/freezing point	Not available.		
Initial boiling point and boiling range	Not available.		
Pour point	Not available.		
Specific gravity	Not available.		
Partition coefficient (n-octanol/water)	Not available.		
Flash point	65.3 °F (18.5 °C) Tag Closed Cup (liquid)		
Evaporation rate	Not available.		
Flammability (solid, gas)	Not applicable.		
Upper/lower flammability or exp	losive limits		
Flammability limit - lower (%)	Not available.		
Flammability limit - upper (%)	Not available.		
Explosive limit - lower (%)	Not available.		
Explosive limit - upper (%)	Not available.		
Vapor pressure	Not available.		
Vapor density	Not available.		
Relative density	Not available.		
Solubility(ies)	Not available.		
Auto-ignition temperature	Not available.		
Decomposition temperature	Not available.		
Viscosity	3.23 cSt @ 68°F (liquid)		
Other information			
Density	0.86 (liquid)		
	10. Stability and Reactivity		
Reactivity	This product may react with strong oxidizing agents.		
Possibility of hazardous reactions	No dangerous reaction known under conditions of normal use.		
Chemical stability	Stable under recommended storage conditions.		
Conditions to avoid	Avoid heat, sparks, open flames and other ignition sources. Avoid temperatures exceeding the flash point. Do not mix with other chemicals.		
Incompatible materials	Strong oxidizing agents.		
Hazardous decomposition products	May include and are not limited to: Oxides of carbon.		
	11. Toxicological Information		
Routes of exposure	Eye, Skin contact, Inhalation, Ingestion.		
Information on likely routes of e			
Ingestion	May cause stomach distress, nausea or vomiting.		
Inhalation	Vapors have a narcotic effect and may cause headache, fatigue, dizziness and nausea. Harmful if inhaled. Prolonged inhalation may be harmful.		
Skin contact	No adverse effects due to skin contact are expected. Non-irritating based on test data.		
Eye contact	Moderate eye irritation		
Symptoms related to the physical, chemical and toxicological characteristics	Symptoms may include stinging, tearing, redness, swelling, and blurred vision. Symptoms of overexposure may be headache, dizziness, tiredness, nausea and vomiting.		
Information on toxicological effe	ects		
Acute toxicity	Harmful if inhaled. Narcotic effects. May cause respiratory irritation.		
·····			

Product	Species	Test Results
IAIC 88494-2 (Wipe) (CAS Mi	xture)	
Acute		
<i>Dermal</i> LD50	Rat	> 5000 mg/kg Toptod
	Ral	> 5000 mg/kg, Tested
Inhalation LC50	Rat	> 2.1 mg/L, Tested
	Nat	> 2.1 mg/L, 163teu
<i>Oral</i> LD50	Rat	> 5000 mg/kg, Tested
components	Species	Test Results
idecyl dimethyl ammonium ch	_	
Acute		
Dermal		
LD50	Rabbit	500 mg/kg
	Rat	> 1000 mg/kg
Inhalation		
LC50	Not available	
Oral		
LD50	Mouse	268 mg/kg
	Rat	329 mg/kg
		84 mg/kg
thanol (CAS 64-17-5)		
Acute		
Dermal		
LD50	Rabbit	> 15800 mg/kg, SIDS initial assessmen report
		report
Inhalation LC50	Cat	85.4 mg/L, 4.5 Hours, ECHA
2030	Cal	-
		43.7 mg/L, 6 Hours, ECHA
	Mouse	> 60000 ppm, 60 Minutes, ECHA
		79.4 mg/L, 134 Minutes, ECHA
	Rat	> 115.9 mg/L, 4 Hours, ECHA
		31623 ppm, 4 Hours, HMIRA
		20000 ppm, 10 Hours, HSDB
		51.3 mg/L, 6 Hours, ECHA
Oral		
LD50	Dog	5.5 g/kg, HSDB
	Guinea pig	5600 mg/kg, HSDB
	Monkey	6000 mg/kg
	Mouse	10500 ml/kg, ECHA
		3450 mg/kg, SAX
	Pig	> 5000 mg/kg, ECHA
	Rat	1187 - 2769 mg/kg, ECHA
		12400 mg/kg, ECHA
		10470 mg/kg, ECHA
		7800 ml/kg, ECHA
kin corrosion/irritation	Prolonged skin contact may cause tempor	ary irritation.
Exposure minutes	Not available.	
Erythema value	Not available.	
Oedema value	Not available.	
Serious eye damage/eye	Causes serious eye irritation.	
rritation		

Corneal opacity value	Not available			
Iris lesion value	Not available	9.		
Conjunctival reddening value	Not available	Э.		
Conjunctival oedema value	Not available	Not available.		
Recover days	Not available	Not available.		
Respiratory or skin sensitizatio	'n			
Respiratory sensitization	Not available	Not available.		
Skin sensitization	This product	is not expected to	cause skin sensitization.	
Mutagenicity	Non-hazardo	ous by WHMIS/OS	SHA criteria.	
Carcinogenicity				
Canada - Manitoba OELs: c	arcinogenicity	/		
ETHANOL (CAS 64-17- IARC Monographs. Overall		Carcinogenicity	Confirmed animal carcin	ogen with unknown relevance to humans.
Ethanol (CAS 64-17-5)			Volume 44, Volume 96, Volume 96, Volume 96, Volume 100	
US - California Proposition		ed date/Carcinog	enic substance	
consumption.	d to be a Prop 6	·		is not an alcoholic beverage for
US. OSHA Specifically Reg Not listed.	ulated Substal	nces (29 CFR 191	0.1001-1050)	
Reproductive toxicity	Not available	2		
Teratogenicity		ous by WHMIS/OS		
•		-		
Specific target organ toxicity - single exposure		Respiratory tract irritation. Narcotic effects.		
Specific target organ toxicity - repeated exposure	Not classifie	d.		
Aspiration hazard	Not available	Not available.		
Chronic effects	Prolonged in	Prolonged inhalation may be harmful.		
		12. Ecologic	al Information	
Ecotoxicity	See below			
Ecotoxicological data				
Components		Species		Test Results
Didecyl dimethyl ammonium chlor	ride (CAS 7173	-51-5)		
Aquatic				
Fish	LC50	White sturgeor transmontanus		0.001 - 0.01 mg/L, 96 hours
Ethanol (CAS 64-17-5)				
Crustacea	EC50	Daphnia		11744.5 mg/L, 48 Hours
Aquatic				
Crustacea	EC50	Water flea (Da	phnia magna)	7.7 - 11.2 mg/L, 48 hours
Fish	LC50			
-				
Persistence and degradability	No data is available on the degradability of this product.			
Bioaccumulative potential	No data available.			
Mobility in soil	No data available. Not available.			
Mobility in general				
Other adverse effects	No other adverse environmental effects (e.g. ozone depletion, photochemical ozone creation potential, endocrine disruption, global warming potential) are expected from this component.			

	13. Disposal Considerations	
Disposal instructions	 Do not contaminate water, food or feed by storage or disposal. PESTICIDE DISPOSAL: Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility. CONTAINER DISPOSAL: Non-refillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash. CONTAINER DISPOSAL: Refillable container. Refill plastic canister with this brand of pesticide product only. Do not reuse this container for any other purpose. Offer plastic container for recycling if available, or discard in trash. Collect and reclaim or dispose in sealed containers at licensed waste disposal site. This material and its container must be disposed of as hazardous waste. Incinerate the material under controlled conditions in an approved incinerator. Do not incinerate sealed containers. Do not allow this material to drain into sewers/water supplies. Do not contaminate ponds, waterways or ditches with chemical or used container. If discarded, this product is considered a RCRA ignitable waste, D001. Dispose of contents/container in accordance with local/regional/national/international regulations. 	
Local disposal regulations	Dispose in accordance with all applicable regulations.	
Hazardous waste code	D001: Waste Flammable material with a flash point <140 F The waste code should be assigned in discussion between the user, the producer and the waste disposal company.	
Waste from residues / unused products	Dispose of in accordance with local regulations. Empty containers or liners may retain some product residues. This material and its container must be disposed of in a safe manner (see: Disposal instructions).	
Contaminated packaging	Empty containers should be taken to an approved waste handling site for recycling or disposal. Since emptied containers may retain product residue, follow label warnings even after container is emptied.	
	14. Transport Information	
Transport of Dangerous Goods (TDG) Proof of Classification	Classification Method: Classified as per Part 2, Sections 2.1 – 2.8 of the Transportation of Dangerous Goods Regulations. If applicable, the technical name and the classification of the product will appear below.	
General	IATA: Passenger and Cargo Aircraft: Limited Quantity: Pkg Inst: Ltd Qty Y341 Passenger and Cargo Aircraft: Limited Quantity: Pkg Inst: Max Net Qty/Pkg 1L	
	Passenger and Cargo Aircraft: Pkg Inst: Pkg Inst 352 Passenger and Cargo Aircraft: Pkg Inst: Max Net Qty/Pkg 1L	
	CARGO Aircraft Only: Pkg Inst 364 CARGO Aircraft Only: Max Net Qty/Pkg 60L	
	IMDG: Special Provisions: 274 Limited Quantity: 1L Excepted Quantity: E2 Packing Instructions: P001 IBC Instructions: IBC02 Portable Tanks and Bulk Containers: Tank Instructions: T7 Provisions: TP2, TP13 Stowage and Handling: Category B, SW2	
U.S. Department of Transportation	ion (DOT)	
Basic shipping requirement		
UN number	UN1992	
Proper shipping name	Flammable liquids, toxic, n.o.s.	
Technical name Hazard class	Ethanol Limited Quantity - US	
Packing group	II	
Special provisions	IB2, T7, TP2, TP13	
Packaging exceptions	<0.3 Gal-Consumer Commodity ORM-D/Limited Quantity	
Packaging non bulk Packaging bulk	202 243	
Transportation of Dangerous G		
Basic shipping requirement		
UN number	UN1992	
Proper shipping name	FLAMMABLE LIQUID, TOXIC, N.O.S.	
Technical name	Ethanol	
Hazard class Packing group	Limited Quantity - Canada II	

Packaging exceptions	16 <1L - Limited Quantity		
IATA/ICAO (Air)			
Basic shipping requiremer	nts:		
UN number	UN1992		
Proper shipping name	Flammable liquid, toxic, n.o.s.		
Technical name Hazard class	Ethanol Limited Quantity - IATA		
Packing group			
ERG code	3HP		
IMDG (Marine Transport)			
Basic shipping requirement			
UN number Broner chinning name	UN1992 FLAMMABLE LIQUID, TOXIC, N.O.S.		
Proper shipping name Technical name	Ethanol		
Hazard class	Limited Quantity - IMDG		
Packing group	II		
Marine pollutant	Yes		
EmS	F-E, S-D		
DOT; IMDG; TDG			
$\langle \rangle$			
ΙΑΤΑ			
$\langle \mathbf{Y} \rangle$			
$\langle \mathbf{Y} \rangle$			
Y			
Y			
Y	15. Regulatory Information		
Canadian federal regulations	15. Regulatory Information This product has been classified in accordance with the hazard criteria of the HPR and the SDS		
Canadian federal regulations	15. Regulatory Information This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR.		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR.		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II).		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood.		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS:		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS: Hazards to humans and Domestic animals.		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS:		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS: Hazards to humans and Domestic animals. CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly		
Canadian federal regulations	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS: Hazards to humans and Domestic animals. CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating and drinking, chewing gum, using tobacco o using toilet. PHYSICAL OR CHEMICAL HAZARDS		
	This product has been classified in accordance with the hazard criteria of the HPR and the SDS contains all the information required by the HPR. This product is not subject to the Hazardous Products Act (HPA) Part II (Hazardous Products) as per paragraph 12(j); Schedule 1 (Non-Application of Part II). This restriction states that Part II does not apply in respect of the sale or importation of anything listed in Schedule 1 which includes any pest control product as defined in subsection 2(1) of the Pest Control Products Act, any explosive as defined in section 2 of the Explosives Act, any cosmetic, device, drug or food, as defined in section 2 of the Food and Drugs Act, any consumer product as defined in section 2 of the Canada Consumer Product Safety Act and any wood or product made of wood. KEEP OUT OF REACH OF CHILDREN. CAUTION PRECAUTIONARY STATEMENTS: Hazards to humans and Domestic animals. CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating and drinking, chewing gum, using tobacco o using toilet.		

Not listed. FreeNues Gases Not isted. FreeUser Control Regulated. Not regulated. Not regulated. Not regulated. US foderal regulated. IS separation Separate Control Regulated. IS a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines. IS for the control Regulated. IS a registered EPA Registration No EPA 88494-2 This chemical is a disrification product registered by the United States Environmental Protection Agency and is subject to certain heading requirements under federal pesticide law. These requirements differ from the classification refere and heazed information required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required on the disrification label is reproduced below. The disrification required information required with scop and water after handing and before eating and diriking, chewing gum, using tobacco or using bolde. FEECALTIONARY STATEMENTS: INTERCA Section 12(b Export NetFication (40 CFR 302.4) Not listed. Superfund Amendments and Resutter Section 2000	Export Control List (CEPA 1	999, Schedule 3)	
Not listed. Not applicable WHMB 2015 Exemptions Not applicable US federal regulations This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines. EPA Registration No EPA 88494-2 EPA Registration No EPA 88494-2 This chemical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to certain labeling requirements diffected repeated by the United States Environmental Protection Agency and is subject to certain labeling requirements diffected refere a pescicie law. These requirements differ from the classification criteria and hazard information required for safety data sheets (SDS), and for workplace labels of non-disinfectant to the hazard information required to the safety data sheets (SDS), and for workplace labels of non-disinfectant to the hazard information required to the safet so ruse. PRECAUTIONC subse moderate evy initiation. Avoid on the Azard information required to the safet so ruse. The Azard to humans and Domestic animals. CAUTIONC causes moderate evy initiation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handing and before eating and diriking, chewing gum, using tobacc or using tobac. Not regulated. Elementation (40 CFR 707, Subpt. D) Not regulated.			
Precursor Control Regulations Not applicable WHMIS 2015 Examptions Not applicable US federal regulations This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines. EPA Registration No EPA 88494-2 EPA Registration No EPA 88494-2 This charmical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to cartin labeling reguirements under federal possicide state. The saver information required on the disinfectant labeling requirements under federal possicide states on on-disinfectant chemicals. The hazard information, including divelsions for used infectant chemicals. The hazard information, including divelsions for used infectant chemicals. The hazard information, including divelsions for user divelsions for tester divelsions for user divelsions for tester di			
Not regulated. Not applicable WHMIS 2015 Exemptions Not applicable US federal regulations This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines. EPA Registration No EPA 88494-2 This chamical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to cortain labeling requirements and harzard information required for safety data sheets (GDS), and for wareplace labels of non-disinfectant chemicals. The horization for use. PRECAUTIONARY STATEMENTS: Harzards to humans and Domestic animals. CAUTION: Cause moderate ey irritation. Avoid contact with eyes or clothing. Wash thoroughly with scap and water after handling and before eating and dimking, chewing gum, using tobacco or using tolia. Mot regulated. ERECAUTIONARY STATEMENTS: Hazard to humans and Domestic animals. CAUTION: Cause moderate ey irritation. Avoid contact with eyes or clothing. Wash thoroughly with scap and water after handling and before eating and dimking, chewing gum, using tobacco or using tolia. Mot regulated. US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not listed. Immediate Harazrd - Yes Superfund Amendments and Reauthorization Act of 1986 (SARA) Hazard oats substance SARA 3102 Extremely No Not regulated. Reactivity Hazard - Yes Pressum Haz		ins	
WHIS 2015 Exemptions Not a regulations US federal regulations This is a regulated EPA product. The product labeling is in compliance with EPA regulations and guidelines. UR Federal regulations EPA Registration No EPA 88494-2 This chemical is a disinfectant product registered by the United States Environmental Protection Apperoy and is subject to certain and hazard information regulators. The hazard information regulator does allowed for advision criteria and hazard information regulators and hazard information regulators and hazard information regulators. The hazard information regulators and hazard information regulators and hazard information regulators. The hazard information regulators and more ables of one-distribution criteria and hazard information regulators and hazard information regulators. The CAUTIONARY STATEMENTS: Hazard to humans and Domostic animals. CAUTION. Causes moderate eye intition. Avoid contact with eyes or clothing. Wash thoroughly with scop and water after handing and before eating and finder exemption gum, using tobacco or using tolet. TSCA Section 12(b) Export Nutification (40 CFR 707, Subpt. D) Not regulator. Not regulator. Examption Act of 1986 (SARA) Superfund Amendments Apperities at Substances (29 CFR 1910.1001-1050) Not listed. Superfund Amendmentare and regulators. Not regulator. SARA 302 Extremely No Pressure Lazard - No Pressure Lazard - No SARA 3113(TRI reporting) Not regulator. Us celliforia H	•		
US federal regulations This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines. EPA Registration No EPA 88494-2. This chemical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to certain labeling requirements and hazard information required for safety data sheets (EUS), and for workplace labels of hom-disinfectant chemicals. The hazard information important information, including directions for use. PRECAUTIONARY STATEMENTS: Hazardo to humans and Domestic animals. CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with scap and water after handling and before eating and drinking, chewing gum, using tobacco or using tolie. TSCA Section 12(b) Export Notification (40 CFR 707, Subpl. D) Not regulated. Superfund Amazmetments and Reauthorization Act of 1986 (SARA) Hazard ox Substance List (40 CFR 707, Subpl. D) Not listed. Immediate Hazard - Yes Superfund Amadments and Reauthorization Act of 1986 (SARA) Hazard ox Substance Hazard categories Immediate Hazard - Yes Pressure Hazard - Yes Delayed Hazard - Yes SARA 313 (TRI reporting) No Not regulated. No US california Hazardo Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Us california Hazard - Yes Listed. SARA 313	-	Not applicable	
Arris chemical is a disinfectant product registered by the United States Environmental Protection Agency and is subject to cartain labeling requirements under federal peation required for safety data sheets (SDS), and for workplace labels of non-disinfectant hemicals. The bazard information required on the disinfectant label is erportude below. The disinfectant label also includes other important information, including directions for use. PRECAUTIONARY STATLEMENTS: Hazards to humans and Domesite animals. CAUTION: Causes moderate use inflation. Avoid contact with eyes or clothing. Wash thoroughly with scap and water after handling and before eating and diriking, chewing gum, using tobacco or using tolat. TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D) Not regulated. TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D) Not regulated. Superfund Amendments and Reaution Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Pressure Hazard - No Pressure Hazard - No Pressure Hazard - No Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Not regulated. SARA 302 Extremely Not regulated. No Reactivity Hazard - No Pressure Hazard - No Pressure Hazard - No Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Pressure Hazard - No Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Pressure Hazard - No Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Pressure Hazard - No Reactivity Hazard - N	-	This is a registered EPA produ	uct. The product labeling is in compliance with EPA regulations and
Agency and is subject to certain labeling 'requirements under federal pecticle law. These requirements differ from the classification criteria and hazard information required on the disinfectant label is oproduced below. The disinfectant label also includes other important information, including directions for use. PRECAUTION/ Causes moderate eye inflation. Avoid contact with eyes or clothing. Wash thoroughly with a subject on the disinfectant label is approduced below. The disinfectant label also includes other important information, including directions for use. PRECAUTION/ Causes moderate eye inflation. Avoid contact with eyes or clothing. Wash thoroughly with a subject was and water after handling and before eating and diriking, chewing gum, using tobacco or using poliet. TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D) Not regulated. CERCLA Hazardous Substance: List (40 CFR 302.4) Not listed. Superfund Amendments and Reauthorization Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No Fire Hazard - Yes Pressure Hazard - No Reactivity Hazard - No Reactivity Nakard - No Reactivity Hazard - No Reactivity Nakard - No Reactivity Nakard - No Reactivity Hazard -		EPA Registration No EPA 88	8494-2
Hazards to humans and Domestic animals. CAUTION: Causes moderate eye initation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating and drinking, chewing gum, using tobacco or using toilet. Flammable. Keep away from heat or open flame. TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D) Not regulated. CERCLA Hazardous Substance List (40 CFR 302.4) Not listed. US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not listed. Superfund Amendments and Re-uthorization Act of 1986 (SARA) Hazard Categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - No Reactivity Hazard - No Reactity Hazard - No Reactit Reactivity Hazard - No Reactivity Hazard -		Agency and is subject to certa requirements differ from the cla sheets (SDS), and for workpla required on the disinfectant lab	in labeling requirements under federal pesticide law. These assification criteria and hazard information required for safety data ce labels of non-disinfectant chemicals. The hazard information bel is reproduced below. The disinfectant label also includes other
TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D) Not regulated. Not regulated. CERCLA Hazardous Substance List (40 CFR 302.4) Not listed. Not listed. Superfund Amendments and Reauthorization Act of 1986 (SARA) Immediate Hazard - Yes Delayed Hazard - Yes Delayed Hazard - No Fire Hazard - No Reactivity Hazard - No SARA 302 Extremely No ressure Hazard - No Reactivity Hazard - No SARA 311/312 Hazardous chemical No chemical No chemical Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. SarA 313 (TRI reporting) Not regulated. Clean Air Act (CAA) Section 112 (Accidental Release Prevention (40 CFR 68.130) Not regulated. No ressure Listed Substance US - Satifornia Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) US - Colifornia Hazardous Substances (Director's): Listed Substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) US - Vousiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Yeas Effects Screening Levels: Listed substance Ethanol (CAS 64-17-5) US - Yeas Effects Screening Levels: Listed substance Ethanol (CAS 64-17-5) US - Yeas Effects Screening Levels: Listed substance Ethanol (CAS 64-17-5) US - Yeas Effects		Hazards to humans and Dome CAUTION: Causes moderate with soap and water after hand	estic animals. eye irritation. Avoid contact with eyes or clothing. Wash thoroughly
Not regulated. CERCLA Hazardous Substance List (40 CFR 302.4) Not listed. US: OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not listed. Superfund Amendments and Reauthorization Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No SARA 302 Extremely No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No SARA 311/312 Hazardous No chemical SARA 311 (TRI reporting) Not regulated. Other federal regulated. Other federal regulated. Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - Ninnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - Ninnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - Ninnesota Haz Subst. Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - Texas Effects Screening Levels: Listed substance Ethanol (CAS 64-17-5) Listed. Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.		Flammable. Keep away from h	neat or open flame.
CERCLA Hazardous Substance List (40 CFR 302.4) Not listed. SUS- OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not listed. Superfund Amendments and Reauthorization Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - No Reactivity Hazard - No Reactiv	TSCA Section 12(b) Export N	lotification (40 CFR 707, Subp	ot. D)
US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050) Not listed. Superfund Amendments and Reactivitation Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - No Reactivity Hazard - No Reactivity Hazard - No Reactivity Hazard - No SARA 302 Extremely hazardous substance No SARA 311/312 Hazardous chemical No SARA 313 (TRI reporting) Not regulated. No Other federal regulations Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Other federal regulations Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. US state regulations Listed. US state regulations Listed. US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Usiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Minnesota Haz Subs: Listed substance Listed. US - Notisiana Spill Reporting: Listed substance Listed. <th>0</th> <th>nce List (40 CFR 302.4)</th> <th></th>	0	nce List (40 CFR 302.4)	
Superfund Amendments and Re-uthorization Act of 1986 (SARA) Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - No Reactivity	US. OSHA Specifically Regu	lated Substances (29 CFR 191	10.1001-1050)
Hazard categories Immediate Hazard - Yes Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No Reactivity Hazard - No Reactivity Hazard - No SARA 302 Extremely hazardous substance No SARA 311/312 Hazardous chemical No SARA 311/312 Hazardous chemical No SARA 313 (TRi reporting) Not regulated. No Citean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. No Other federal regulations Listed Ciean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. No US state regulations Listed US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) List			
Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No Reactivity Hazardows Reactivity Hazardows Reactivi	-	-	RA)
hazardous substance SARA 311/312 Hazardous No chemical SARA 313 (TRI reporting) Not regulated. Other federal regulations Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Clean Air Act (CAA) Section 112 (r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Listed. US - Listed Safety Act: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Hazard categories	Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No	
chemical SARA 313 (TRI reporting) Not regulated. Other federal regulations Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Lis		No	
Not regulated. Other federal regulations Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) Listed substance Ethanol (CAS 64-		No	
Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List Not regulated. Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) Listed substance Ethanol (CAS 64-17			
Not regulated. Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed substance	Other federal regulations		
Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130) Not regulated. US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Clean Air Act (CAA) Section	112 Hazardous Air Pollutants	(HAPs) List
US state regulations US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed substance Ethanol (CAS 64-17-5) Listed substance Ethanol (CAS 64-17-5) Listed substance	8	112(r) Accidental Release Pre	evention (40 CFR 68.130)
US - California Hazardous Substances (Director's): Listed substance Ethanol (CAS 64-17-5) Listed. US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Not regulated.		
Ethanol (CAS 64-17-5)Listed.US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5)Listed.US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5)Listed.US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5)Listed.US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5)Listed.US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5)Listed.US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5)Listed.	US state regulations		
US - Illinois Chemical Safety Act: Listed substance Ethanol (CAS 64-17-5) US - Louisiana Spill Reporting: Listed substance Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.		ubstances (Director's): Listed	substance
US - Louisiana Spill Reporting: Listed substance Listed. Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Listed. Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) Listed. US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed. Listed.	US - Illinois Chemical Safety	Act: Listed substance	Listed.
Ethanol (CAS 64-17-5) Listed. US - Minnesota Haz Subs: Listed substance Listed. Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed. Listed.		and listed substance	
Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Ethanol (CAS 64-17-5)		Listed.
Ethanol (CAS 64-17-5) US - Texas Effects Screening Levels: Listed substance Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Ethanol (CAS 64-17-5)		Listed.
Didecyl dimethyl ammonium chloride (CAS 7173-51-5) Listed.	Ethanol (CAS 64-17-5)		
		-	Listed
	Ethanol (CAS 64-17-5)		Listed.


US. Massachusetts RTK - S	ibstance List	
Ethanol (CAS 64-17-5)		
US. New Jersey Worker and	Community Right-to-Know Act	
Not regulated.		
US. Pennsylvania Worker ar	d Community Right-to-Know Law	
Ethanol (CAS 64-17-5)		
US. Rhode Island RTK		
Ethanol (CAS 64-17-5)		
US. California Proposition 6	5	
Ethanol is not considered consumption.	to be a Prop 65 component for this proc	luct. The product is not an alcoholic beverage for
US - California Proposit	ion 65 - CRT: Listed date/Carcinogeni	c substance
Ethanol (CAS 64-17-	5) Listed:	April 29, 2011
Υ.		July 1, 1988
US - California Proposit	ion 65 - CRT: Listed date/Developmer	ital toxin
Ethanol (CAS 64-17-	5) Listed:	October 1, 1987
ventory status		
Country(s) or region	Inventory name	On inventory (yes/no)
Canada	Domestic Substances List (DSL)	Ye
Canada	Non-Domestic Substances List (NDSL) No
United States & Puerto Rico	Toxic Substances Control Act (TSCA)	Inventory Ye

*A "Yes" indicates that all components of this product comply with the inventory requirements administered by the governing country(s)

LEGEND	
Severe	4
Serious	3
Moderate	2
Slight	1
Minimal	0

Disclaimer

16	. Othe	r Information	



The information in the sheet was written based on the best knowledge and experience currently available. Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since conditions of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirements of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document.

	in this document.	
Issue date	18-April-2019	
Version #	03 revision of version issues on 18-February-2019	
Effective date	18-April-2019	
Prepared by	North American Infection Control Ltd.	
Other information	For an updated SDS, please contact the supplier/manufacturer listed on the first page document.	

of the