

Surface Disinfectant Wipes and Spray

A fast, effective and easy way to clean and disinfect





An Effective Solution for Complete Disinfection

Kills a broad spectrum of pathogens in just 1 minute

- Efficacy tests have demonstrated that Monarch Surface Disinfectants are an effective, virucide, bactericide, tuberculocide and fungicide in the presence of organic soil.
- One-step, cleaner & disinfectant saves time.

Fast and efficient application

- One minute contact time allows for faster patient transition.
- Compatible with virtually all surfaces in the dental office.
- Helps prevent cross-contamination from contact with treated surfaces.

Gentle

- EPA approved with a low toxicity rating.
- Non-irritating to skin, Non-corrosive and non-staining.
- Fragrance-free and bleach-free formula doesn't leave an offensive odor and dries without a residue.
- Safe to use on dental chair upholstery made of non-porous vinyl.



Because Monarch products are easy to use and produce predictable results, any of our team members can step in and help turn over the operatory at any time.

Carolyn Miskell, Practice Director

A study performed by the **Dental Advisor** examined the effects caused by an accelerated 5-year interaction between *Ultraleather* and *Dauphine* dental chair upholstery and *Monarch Surface Disinfectant Wipes*. No visible changes in appearance were observed for either fabric treated with *Monarch Surface Disinfectant Wipes*.





Figure 1a and b: Microscopic view of a) treated and b) untreated Ultraleather upholstery Source: Dental Advisor Research Report Number 97, October 2016






Figure 2a and b: Microscopic view of a) treated and b) untreated Dauphine upholstery

Order Information

Order No.	Product Description
H6186	Large Monarch Surface Disinfectant Wipes / 6"x6.75"; 160 ct. per canister
H6171	X-Large Monarch Surface Disinfectant Wipes / 7"x9"; 100 ct. per canister
H6110-N	Monarch Surface Disinfectant Spray / 32 oz. Spray Bottle
H6182	Monarch Dry Wipes Intro Kit / 2 canisters with dry wipes, plus 4 refills (100 ea. 7"x9")
H6183	Monarch Dry Wipes Refill Kit / 12 dry wipe refills (100 ea. 7"x9")



CORPORATE HEADQUARTERS 1295 Walt Whitman Road Melville, New York 11747- 3062 800-247-8324



Monarch is a registered trademark of Air Techniques, Inc. © Copyright 2022 PN 9922-618 RevB

