



Small component. Big Effect.
Switch to hygienic compressed air now.



**AIR
TECHNIQUES** equipped for life


Small particles – big risk.

Bacteria, fungi and viruses are not just invisible to the human eye – many microorganisms are too small for dental compressed air filter systems to deal with. You can rely on the Air Techniques Virus-Bacteria filter to protect you against the risk of contamination. A retention rate of **99.994%*** was confirmed with **test virus PhiX174**.

Regardless of particle size:
Retention rate of

99.99995%*

This corresponds to filter class
ULPA U16 according to EN 1822-1:2019-10 /
ISO 65U according to ISO 29463-3:2018



*Test results from 29.06.2020 recorded by the SGS Institut Fresenius GmbH.

Contaminated compressed air is history!

With the Virus-Bacteria Filter from Air Techniques.



- ✓ Switch to an Air Techniques Virus-Bacteria Filter and benefit from hygienic compressed air
- ✓ The Virus-Bacteria Filter can be retrofitted on all AirStar and AirStar NEO compressors



- ✓ Included as standard on all Air Techniques AirStar and AirStar NEO compressors since May 2020

- ✓ AirStar and AirStar Neo preventative maintenance kits



CORPORATE HEADQUARTERS
1295 Walt Whitman Road
Melville, New York 11747- 3062
800-247-8324



AIR
TECHNIQUES equipped for life