



# DIFFERENT AT FIRST SIGHT

FLOWSTAR SCAVENGING CIRCUITS / FLOWSTAR DOUBLE NASAL HOODS

## FLOWSTAR

Bringing Nitrous Oxide / Oxygen Sedation  
to the next level with multiple innovations

**AIR**  
**TECHNIQUES** equipped for life

# Nitrous Oxide / Oxygen Sedation

From Air Techniques – The Market Leader in Air & Vac



Air Techniques is the pioneer in air compressors and vacuum systems since 1962. We are grateful that dentists trust in our products for almost 60 years. This trust made us the most sold brand for utility room products in the USA and in Canada.

FlowSTAR Nitrous Oxide / Oxygen scavenging systems benefit from a high quality vacuum system. With the help of a reliable vacuum system it can be administered safely within your practice. FlowSTAR Nitrous Oxide / Oxygen Sedation is the next natural fit for Air Techniques.

## Different at First Sight




Adult



Pediatric

# Innovations at a Glance



## Patient Comfort



- Soft inner hood improves seal and patient comfort
- Clear outer hood sits closer to the patient's face, increasing scavenging efficiency of exhaled gases and allowing visual monitoring of patient's breathing
- Latex-free material protects patients and staff who have latex allergies


## Safety for Staff



- Double hood design minimizes ambient nitrous oxide in the operator
- FLOWSTAR Double Nasal Hoods utilize an oval exhalation valve that directs exhaled gases to the built-in scavenger "flow channels", maximizing efficiency US patented\*

## Compatible


- Universal fit with other double tubing scavenging circuits allows for use of your existing systems with the new Air Techniques FLOWSTAR Double Nasal Hoods



## Size Matters



- **Size Matters!** Lower profile provides improved "sight line" visibility for the dental staff and therefore allows better access to the oral cavity
- **Accurate fit!** The softer inner hood seals to the face to ensure proper gas delivery to the patient



Lower profile FlowStar Double Nasal Hood as compared to other available double nasal hoods

Competitor 1  
Competitor 2

# FLOWSTAR Double Nasal Hoods

- Colorful and lightly scented single use hoods engage and delight patients of all ages
- Pediatric and adult sizes available in a variety of scents

12 hoods individually wrapped in a box

Part Number	Product description and contents	Scents
NHA12BD	FLOWSTAR Double Nasal Hood Single Use, Adult, Blueberry Dance Scent	 Blueberry Dance
NHP12BD	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Blueberry Dance Scent	 Blueberry Dance
NHA12SD	FLOWSTAR Double Nasal Hood Single Use, Adult, Strawberry Delight Scent	 Strawberry Delight
NHP12SD	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Strawberry Delight Scent	 Strawberry Delight
NHP12BB	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Bouncy Bubble Scent	 BOUNCY BUBBLE
NHA12CV	FLOWSTAR Double Nasal Hood Single Use, Adult, Chilla Vanilla Scent	 Chilla Vanilla
NHP12CV	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Chilla Vanilla Scent	 Chilla Vanilla
NHA12PG	FLOWSTAR Double Nasal Hood Single Use, Adult, Pure Ghost Non Scented	 PURE GHOST
NHP12PG	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Pure Ghost Non Scented	 PURE GHOST

Compatible with Accutron, Porter and Belmed double tubing scavenger circuits



12 hoods individually wrapped / 12 CO2 adapters in a box

Part Number	Product description and contents
NHA12BD-CO2	FLOWSTAR Double Nasal Hood Single Use, Adult, Blueberry Dance Scent, 12 pack of CO2 adapters
NHP12BD-CO2	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Blueberry Dance Scent, 12 pack of CO2 adapters
NHA12SD-CO2	FLOWSTAR Double Nasal Hood Single Use, Adult, Strawberry Delight Scent, 12 pack of CO2 adapters
NHP12SD-CO2	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Strawberry Delight Scent, 12 pack of CO2 adapters
NHP12BB-CO2	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Bouncy Bubble Scent, 12 pack of CO2 adapters
NHA12CV-CO2	FLOWSTAR Double Nasal Hood Single Use, Adult, Chilla Vanilla Scent, 12 pack of CO2 adapters
NHP12CV-CO2	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Chilla Vanilla Scent, 12 pack of CO2 adapters
NHA12PG-CO2	FLOWSTAR Double Nasal Hood Single Use, Adult, Pure Ghost Non Scented, 12 pack of CO2 adapters
NHP12PG-CO2	FLOWSTAR Double Nasal Hood Single Use, Pediatric, Pure Ghost Non Scented, 12 pack of CO2 adapters

airtechniques.com

# Nitrous Oxide / Oxygen Scavenging Circuits

- Fully autoclavable scavenging circuit provides for safer delivery and monitoring

FLOWSTAR double hood scavenging circuit with single use nasal hoods starter kit

Part Number	Product description and contents
NSC-001	FLOWSTAR Double Hood Scavenging Circuit including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight) without vacuum control and vacuum tubing
NSC-VCB	FLOWSTAR Double Hood Scavenging Circuit w/ Vacuum Control Block including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight)

Compatible with most nitrous oxide / oxygen flowmeters



## Autoclavable Hoods

Part Number	Product description and contents
NHAA101	Adult, FLOWSTAR Double Nasal Hood Autoclavable Clear, single liner, non scented
NHAP101	Pediatric, FLOWSTAR Double Nasal Hood Autoclavable Clear, single liner, non scented



DISTRIBUTED BY  
AIR TECHNIQUES  
1295 Walt Whitman Road  
Melville, New York 11747- 3062  
800-247-8324

© Copyright 2024  
PN 9922-600 Rev. J

**AIR TECHNIQUES** equipped for life