

Scavenging Circuits / Double Nasal Hoods

FLOWSTAR

Bringing Nitrous Oxide / Oxygen Sedation to the next level with multiple innovations



Nitrous Oxide / Oxygen Sedation

From Air Techniques - The Market Leader in Air & Vac



Air Techniques is the pioneer in air compressors and vacuum systems since 1962. We are grateful that dentists trust in our products for almost 60 years. This trust made us the most sold brand for utility room products in the USA and in Canada.

FlowStar Nitrous Oxide / Oxygen scavenging systems benefit from a high quality vacuum system. With the help of a reliable vacuum system it can be administered safely within your practice. FlowStar Nitrous Oxide / Oxygen Sedation is the next natural fit for Air Techniques.


Different at First Sight





Adult Pediatric

Innovations at a Glance



Safety for Staff


- Double hood design minimizes ambient nitrous oxide in the operatory
- Air Techniques Nasal Hoods utilize an oval exhalation valve that directs exhaled gases to the built-in scavenger "flow channels", maximizing efficiency US patented

Patient Comfort

- Soft inner hood improves seal and patient comfort
- Clear outer hood sits closer to the patient's face, increasing scavenging efficiency of exhaled gases and allowing visual monitoring of patient's breathing
- Latex-free material protects patients and staff who have latex allergies

Compatible


 Universal fit with other double tubing scavenging circuits allows for use of your existing systems with the new Air Techniques Double Nasal Hoods





Size Matters

- Size Matters! Lower profile provides improved "sight line" visibility for the dental staff and therefore allows better access to the oral cavity
- Accurate fit! The softer inner hood seals to the face to ensure proper gas delivery to the patient



Air Techniques lower profile Double Nasal Hood compared to other available double nasal hoods

··Competitor 1 ··Competitor 2

Double Nasal Hoods

- Colorful and lightly scented single use hoods engage and delight patients of all ages
- Pediatric and adult sizes available in a variety of scents

12 hoods individually wrapped in a box

Part Number	Product description and contents	Scents
NHA12BD	Double Nasal Hood Single Use, Adult, Blueberry Dance Scent	Blueberry Dance
NHP12BD	Double Nasal Hood Single Use, Pediatric, Blueberry Dance Scent	Surevery Sunce
NHA12SD	Double Nasal Hood Single Use, Adult, Strawberry Delight Scent	Standard Dalet
NHP12SD	Double Nasal Hood Single Use, Pediatric, Strawberry Delight Scent	- Strawberry Delight
NHP12BB	Double Nasal Hood Single Use, Pediatric, Bouncy Bubble Scent	BOUN(Y BUBBLE
NHA12CV	Double Nasal Hood Single Use, Adult, Chilla Vanilla Scent	Ghilla Vanilla
NHP12CV	Double Nasal Hood Single Use, Pediatric, Chilla Vanilla Scent	
NHA12PG	Double Nasal Hood Single Use, Adult, Pure Ghost Non Scented	PURE GROST
NHP12PG	Double Nasal Hood Single Use, Pediatric, Pure Ghost Non Scented	July Mark Siles

12 hoods individually wrapped / 12 CO2 adapters in a box

12 110 0 0 0 11 0 11 0 0 0 11 0 11 0 11		
Part Number	Product description and contents	
NHA12BD-CO2	-CO2 Double Nasal Hood Single Use, Adult, Blueberry Dance Scent, 12 pack of CO2 adapters	
NHP12BD-CO2	CO2 Double Nasal Hood Single Use, Pediatric, Blueberry Dance Scent, 12 pack of CO2 adapters	
NHA12SD-CO2	Double Nasal Hood Single Use, Adult, Strawberry Delight Scent, 12 pack of CO2 adapters	
NHP12SD-CO2	Double Nasal Hood Single Use, Pediatric, Strawberry Delight Scent, 12 pack of CO2 adapters	
NHP12BB-CO2	Double Nasal Hood Single Use, Pediatric, Bouncy Bubble Scent, 12 pack of CO2 adapters	
NHA12CV-CO2	Double Nasal Hood Single Use, Adult, Chilla Vanilla Scent, 12 pack of CO2 adapters	
NHP12CV-CO2	Double Nasal Hood Single Use, Pediatric, Chilla Vanilla Scent, 12 pack of CO2 adapters	
NHA12PG-CO2	Double Nasal Hood Single Use, Adult, Pure Ghost Non Scented, 12 pack of CO2 adapters	
NHP12PG-CO2	Double Nasal Hood Single Use, Pediatric, Pure Ghost Non Scented, 12 pack of CO2 adapters	



CO2 adapter

airtechniques.com

Nitrous Oxide / Oxygen Scavenging Circuits

• Fully autoclavable scavenging circuit provides for safer delivery and monitoring

Double hood scavenging circuit with single use nasal hoods starter kit

Part Number	mber Product description and contents	
NSC-001	Double Hood Scavenging Circuit including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight) without vacuum control and vacuum tubing	
NSC-VCB	Double Hood Scavenging Circuit w/ Vacuum Control Block including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight)	

Autoclavable Hoods

Part Numb	Product description and contents	
NHAA101	Adult, Double Nasal Hood Autoclavable Clear, single liner, non scented	
NHAP101 Pediatric, Double Nasal Hood Autoclavable Clear, single lir non scented		

