


*Thank You and Congratulations* on your purchase of the Monarch Safety Tray The tray, made of resistant polypropylene, provides a convenient system for cleaning instruments in the dental clinic, practice and laboratory settings.

Features Include:

- Safety Tray made of resistant polypropylene.
- Filling indicators 1, 2 and 3 quarts (liters).
- Ergonomically shaped carrying handles to enable box to be carried safely.
- Tilt resistant draining position for the sieve insert.
- Specially shaped slits allow complete draining and prevent the instruments from becoming entangled.
- Easy care because of machine washable materials.



## Directions for Use:



www.airtechniques.com 1295 Walt Whitman Road | Melville, New York 11747- 3062 Phone: 800-247-8324 | Fax: 888-247-8481