

Reliably Efficient



EASY-TO-USE LARGE TOUCHSCREEN



TWO SLOTS FOR EFFICIENT WORKFLOW



CONNECTS VIA WI-FI, LAN OR USB

WORK SMARTER


PREVIEW IMMEDIATELY



Check X-ray results quickly on the ScanX Duo Touch's high resolution screen, so you can be confident you have a quality scan before finalizing.

DIGITAL IMAGING

The Air Techniques VisionX Imaging Software, now powered with AI, automatically rotates, optimizes and saves all images. This combined with our new IDX PSPs brings highest resolution currently possible. With automatic adjustable contrast, sharpness and other possible diagnostic filter settings, all make it easier to diagnose and to explain potential treatment plans to patients.



SIMULTANEOUS SCANNING



With its flexible connectivity options, scan jobs can be sent to the ScanX Duo Touch from any workstation. It can also handle multiple scan jobs simultaneously requested from different stations.

HI-RES IMAGING FOR ALL INTRAORAL FORMAT IDX PSPS



Size 0




Size 1



Size 2



Size 3



Size 4

ScanX Duo Touch Specifications

Screen	7" touchscreen, 800 x 480 pixels, 16.7 million colors
Plate Sizes	All intraoral formats, sizes 0, 1, 2, 3, 4 IDX PSPs
Theoretical Resolution Lp/mm	40, 2000 dpi
Greyscale (bits)	65,536 (16)
Weight	17 lbs
Dimensions W x D x H	9.2" x 11.6" x 14.9" / 25.0 x 28.5 x 39.1cm
Interfaces	WiFi, LAN, USB

