

UTILITY ROOM ACCESSORIES

# Utility Room Accessories

Take your equipment's performance to the next level.





## Increase The Performance Of Your Equipment

Accessories for your utility room equipment are just as important as the equipment itself. Accessories offer a wide range of benefits such as: space savings, environmental responsibility, energy savings and more. Along with proper cleaning and maintenance, accessories can help keep your equipment running smoothly and extend the life of your gear.

Protect your office and investments with Air Techniques accessories. Office downtime from equipment servicing is reduced when preventative accessories are used. Keeping your office functioning smoothly means reducing down-time so you are able to see more patients.

Create a more efficient office with accessories that increase performance of your compressor or vacuum and utilize only the power you need. With just a few simple product additions, you can ensure that you are getting the highest quality performance from the products you rely on to keep your practice running.

### **Evacuation System Cleaner**

Monarch CleanStream cleans vacuum systems and prevents debris buildup.

- Prolonged vacuum life and uninterrupted workflow.
- Easy to use dispenser cleans two operatories with one bottle.
- Works for both wet and dry vacuum systems.
- Non-foaming, enzyme based, neutral pH.
- Environmentally friendly.





**Amalgam Separator** 

Acadia Plus safely and effectively collects amalgam waste.

- Supports up to 10 operatories
- Easy, no-hassle installation
- Free-standing or Wall-mounted
- Translucent container provides a clear view into filter
- Filter twists on/off for less-than-a-minute replacements
- Prepaid shipping and packaging with the Acadia Plus Recycling Kit



Acadia Plus Clear Filter



- Installed in easily accessible locations.
- Ability to completely shut off water supply.
- 24 volts or 115 volt models offered.
- Available with or without filter.



### Vacuum Equalizer

Monitors vacuum demand throughout the day to provide optimal suction for your practice.

- Automatically switches pumps on and off as needed.
- Reduces energy use and extends service life.
- Boosts vacuum pump efficiency for up to 4 pumps.
- Controls all your pumps with one remote panel switch.

### Hydromiser

Air/Water separation and water recycling system for wet-ring vacuums.

- Properly vents noxious gases like nitrous from your office.
- Recycles up to 75% of water used.
- Saves over 65,000 gallons of water per year.
- Operates without any electricity or moving parts.











### Vision Monitor

Intelligent utility room cloud based monitoring system.

- IoT Condition Monitoring system with web-based user interface (on web and mobile App: Apple & Android)
- Alarms and (pre-) warnings are reported immediately to the desktop and by email
- Remote access to practice equipment via Cloud: easy, safe and secure device configuration



Connect Box

Connect your AirStar Classic compressor to the Vision Monitor IoT platform to ensure you're always up and running. By using the Connect Box you can add any of the

AirStar Classic models, or other Air Techniques equpiment, into Vision Monitor's cloud monitoring system to ensure your entire utility room gets connected.



Control your entire utility room using remote switches direct from your office.

- Installed centrally located in your office.
- Customizable to fit your practice needs.
- Streamlined modern design with a stainless steel cover-plate.
- All four control panels fit over a standard 2 gang electrical box.

# Ur office.



### Remote Air Intake

Delivers fresh air to the compressor heads.

- Enables you to comply with local building codes.
- Reduces compressor head noise by 2-3 decibels.
- Can be mounted vertically or horizontally.
- Helps keep compressor motors cool in restrictive environments.

### Utility Stack Racks

Optimize space and organize your utility room.

- Available in two sizes: Mini and XL.
- 3 height adjustments for upper shelf.
- Adjustable top rails to accommodate all equipment sizes.
- Accommodates either air atop vacuum system or vice versa.





#### **Buck Boost Transformer**

Regulates the line voltage to your equipment.

- If voltage is low, it will boost to a higher voltage by 10%.
- If voltage is high, it will buck to a lower voltage by 10%.

## Specifications

Product	Part Number	Description
Acadia Plus	APLUS	Acadia Plus Amalgam Separator
	APREC	Acadia Plus Collection Container Cartridge & Recycle Kit
	APREP	Acadia Plus Collection Container Cartridge Replacement
	APKIT	Acadia Plus Collection Container Recycle Kit
Stack Racks	56600	Utility Stack Rack Mini
	56650	Utility Stack Rack XL
	56635	Drip Tray Kit for Utility Stack Rack Mini
	56655	Drip Tray Kit for Utility Stack Rack XL
	54145	Stack Rack Wall Mount Anchor Bracket Kit
	56645	VacStar 20 & 40 Top Shelf Mount Kit
Monarch	57850	Monarch CleanStream 2.5 Liter Bottle
	57665	Monarch CleanStream Dispenser System
	57660	Monarch CleanStream Starter Kit
Remote Control Panels	53111	Remote Control Panel – 1 Switch
	53251	Remote Control Panel – 2 Switch
	53250	Remote Control Panel – 3 Switch
	53133	Remote Control Panel – 4 Switch
Buck/Boost Transformer	67006	Buck/Boost Transformer 120 VAC - 2.5 kVA
	67002	Buck/Boost Transformer 220 VAC - 3.9 kVA
	67000-1	Buck/Boost Transformer 220 VAC - 7.8 kVA
VacStar Accessories	H-2	Hydromiser for VacStar 20
	H-4	Hydromiser for VacStar 40
	56041	Hydromiser for VacStar 50
	56042	Hydromiser for VacStar 80
	55540	Air/Water Separator for VacStar
	56200	Vacuum Equalizer
	56230	Remote Switch Panel for Vacuum Equalizer
Remote Water Control	53020	Remote Water Control Valve – ¾" Pipe / 24V - w/ Filter & Transformer
	53020-1	Remote Water Control Valve – ¾" Pipe / 115V - w/ Filter
	53021	Remote Water Control Valve – ¾" Pipe / 24V w/ Transformer
	53021-1	Remote Water Control Valve – ¾" Pipe / 115V
	53170	Remote Water Control Valve – 1" Pipe / 24V - w/ Filter & Transformer
	53171	Remote Water Control Valve – 1" Pipe / 115V - w/ Filter
	53173	Remote Water Control Valve – 1" Pipe / 24V - w/ Transformer
	53174	Remote Water Control Valve – 1" Pipe / 115V
AirStar Accessories	85491	Remote Air Intake – AS10/NEO, AS12
	85492	Remote Air Intake – AS21/NEO, AS22/NEO, AS30/NEO
	87361	Remote Air Intake – AS40/NEO
	85493	Remote Air Intake – AS50/NEO
	85494	Remote Air Intake – AS70/NEO
	85511	Remote Air Intake – AS100
	89831	Air Intake Filter Replacement (1 p/cylinder)
	85694	AirStar Tandem Install Kit (AirStar Classic only)
	85690	AirStar Multiple Compressor Timer Relay Kit with Transformer (up to 4 AirStars)



airtechniques.com

### 60 YEARS OF PROVEN PERFORMANCE

60 years ago, Air Techniques introduced its first dental air compressor. Today, Air Techniques air compressors and vacuum systems are the standard for dental utility equipment. They are used to power more dental offices in North America than all other air compressors and vacuum systems combined.

Air Techniques is the leader in dental air compressor and vacuum systems. Nearly 7 of 10 air compressors and 6 of 10 vacuum systems sold today are made by Air Techniques.









