

Dental Unit Waterline Cleaner

Proven formula for the cleaning and control of microbial contamination in dental unit waterlines



Chlorhexidine based Shock and Maintenance Treatment

Exceeds EPA regulatory standards for drinking water

- Reduces and maintains CFU (colony forming units) counts of <200 CFU/mL with regular use in dental unit waterlines (DUWL).
- Removes biofilm from DUWL containing odor causing bacteria.

Non-corrosive and gentle to equipment

- pH-neutral formula provides initial shock and weekly maintenance treatment in one solution without the use of harsh chemicals.
- Fresh, minty taste.

No mixing or diluting required

• Ready-to-use solution can be used with any type of independent water bottle system. Doesn't leave a residue behind.

Fast, effective, and efficient application

- Purge DUWL, add Lines Cleaner, flush next day.
- Does not affect bond strength.

Every dental practice should be using Monarch Lines Cleaner," which is "an extremely economical solution, especially considering the importance of safety for our patients.





Before

Build-up on internal

walls in just 3 weeks.

After

Clean, smooth walls.

No build-up.

Edward Lew, DMD

Order Information

Order No. Product Description

H6340 Monarch Lines Cleaner 16.9 oz. / 500 mL Bottle; 10 bottles per shipping case



