## AIRSTAR & AIRSTAR NEO

## Semi-Annual/Annual Filter Change and Preventive Maintenance

- Based on Air Quality, Dust and (if using) Remote Air Intake (Every 6-12 Months)
  - Dirty Air Filters can increase Tank Charge-up Times by 150% - 250%



## PART# PREVENTATIVE MAINTENANCE KITS

87351	<b>Preventative Maintenance Kit for AS10 / AS10 NEO &amp; AS12</b> Includes: Qty. 1 compressor head intake filter, Qty. 1 Virus-Bacteria filter 0.01 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 maintenance record card.	
87352	Preventative Maintenance Kit for AS21/AS21 NEO, AS22/AS22 NEO & AS30/AS30 NEO Includes: Qty. 2 compressor head intake filter, Qty. 1 Virus-Bacteria filter 0.01 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 filter bowl, Qty. 1 maintenance record card.	
87355	<b>Preventative Maintenance Kit for AS40 &amp; AS40 NEO</b> Includes: Qty. 3 compressor head intake filter, Qty. 1 Virus-Bacteria filter 0.01 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 filter bowl, Qty. 1 maintenance record card.	
87353	<b>Preventative Maintenance Kit for AS50 &amp; AS50 NEO</b> Includes: Qty. 4 compressor head intake filter, Qty. 1 Virus-Bacteria filter 0.01 microns, Qty. 1 membrane bottom filter 35 microns, Qty. 1 top cover o-ring, Qty. 1 bottom filter bowl o-ring, Qty. 1 filter bowl, Qty. 1 maintenance record card.	
87354	<b>Preventative Maintenance Kit for AS70 &amp; AS70 NEO</b> Includes: Qty. 6 compressor head intake filter, Qty. 2 Virus-Bacteria filter 0.01 microns, Qty. 2 membrane bottom filter 35 microns, Qty. 2 top cover o-ring, Qty. 2 bottom filter bowl o-ring, Qty. 2 filter bowl, Qty. 2 maintenance record card.	
87356	<b>Preventative Maintenance Kit for AS100</b> Includes: Qty. 8 compressor head intake filter, Qty. 2 Virus-Bacteria filter 0.01 microns, Qty. 2 membrane bottom filter 35 microns, Qty. 2 top cover o-ring, Qty. 2 bottom filter bowl o-ring, Qty. 2 filter bowl, Qty. 2 maintenance record card.	

