CATEGORY SPOTLIGHT



FREE SAMPLE

Scan the QR code below to fill out a form and receive a free sample of Monarch CleanStream Evacuation System Cleaner.



AIR TECHNIQUES

Monarch Lines Cleaner & Monarch CleanStream

aintaining safety and efficiency within a dental practice is non-negotiable. However, dental unit waterlines and evacuation systems pose a dual threat. Untreated waterlines jeopardize the water quality essential for various dental procedures, potentially endangering patient safety. Meanwhile, the buildup of debris within evacuation lines can hinder suction performance while increasing the risk of cross-contamination.

Fortunately, Monarch Lines Cleaner and CleanStream Evacuation System Cleaner by Air Techniques ensure dental practices need not compromise on either front. By effectively treating waterlines and evacuation lines, these solutions extend the life of equipment and help facilitate safer, more efficient dental visits.

Monarch

Proven Formula for Waterline Treatment

The Monarch Lines

Cleaner is a noncorrosive, ready-to-use solution that provides initial shock and weekly maintenance treatment in 1 product. Its fast, effective, and efficient application replaces daily water cleaning tablets, leaving no residue in the lines, and removes odorcausing bacteria and biofilm

from waterlines without harsh chemicals.

The chlorhexidine-based formula is also pH neutral and gentle on office equipment.

Compatible with any independent water bottle system, the cleaner has a fresh, minty taste

and does not affect bond strength. One 500mL bottle of Monarch Lines Cleaner can treat a single operatory for several weeks, making it an economical waterline treatment option, too.



Optimal Evacuation System Performance

Designed for daily cleaning, Monarch CleanStream Evacuation System Cleaner removes vacuum fluids and debris, preventing future buildup and improving suction. The concentrated, nonfoaming formula cleans and deodorizes evacuation lines for all wet and dry vacuum systems, and an 84.5-oz bottle yields 125 cleanings. The solution lasts up to 5 days, if premixed.

Monarch CleanStream's innovative dispenser system with Vortex technology atomizes the mixture to clean all inner tubing walls, affording easy, no-mess mixing and cleaning of 2 operatories. The dispenser automatically halts suction once the first cleaning process is complete. The closed dispensing system ensures spill-free transportation to the next unit.

Monarch CleanStream is noncorrosive, pH neutral, and fully compliant with Best Management Practices for amalgam waste.

FOR MORE INFORMATION:

800.247.8324 airtechniques. com