

Film Processor Maintenance

IMAGING | UTILITY ROOM | MERCHANDISE

YOU DEPEND ON RADIOGRAPHIC FILMS FOR DIAGNOSIS AND OFTEN EMPLOY THEM IN THE PRESENTATION OF TREATMENT PLANS OR AS A PATIENT EDUCATION AID.

YOU REALIZE THE IMPORTANCE OF CONSISTENT, ARTIFACT FREE FILM PROCESSING.

- Regular maintenance of your processor with the right maintenance products ensures best quality radiographs and longest service life.
- Air Techniques' maintenance products are specially formulated for maximum cleaning efficiency and ease-of-use.
- Build-up residues are easily removed and staining is reduced, protecting films from roller deposits.
- All Air Techniques' maintenance products are environmentally safe; they do not contain any heavy metals.

FORMULA 2000 PLUS

The non-toxic and non-hazardous cleaning concentrate effectively removes deposits and stains from roller transports, processors and dip tanks.



Film Processor Maintenance

Maintenance for Roller Transport Processors with Formula 2000 PLUS and Spray 2000 PLUS

Quarterly Cleaning with Formula 2000 PLUS

- Reduces staining and prolongs the life of your processor.
- Non-toxic and non-hazardous.
- Effectively dissolves deposits, dirt, algae and residue.
- Deep cleaning of the Roller Transport Processors takes only 15 minutes.

Weekly Cleaning with Spray 2000 PLUS

- Removes residue and protect films from roller deposits.
- Available in a convenient, ready to use spray bottle.
- No mixing necessary - just spray, sponge and rinse under running water.
- For use in A/T 2000 XR.



(Left): Formula 2000 Plus Cleanser Twin Pack, PN 43945
(Right): Spray 2000 Plus, PN 43965

Maintenance for Peri-Pro Transports with Peri-Pro Transport Cleaning Kit

Peri-Pro Transport Cleaning Kit with Formula 2000 PLUS

- Reduces staining, decreases normal maintenance time and prolongs the life of your Peri-Pro.
- Easy and convenient way to care for your Peri-Pro.
- Specially designed to hold the transport with just the right amount of Formula 2000 PLUS and water.
- Tub is completely self-contained, portable and requires little counter space.
- Deep cleaning only takes 30 minutes.
- Can be used for all models of Peri-Pro.



Peri-Pro Transport Cleaning Kit, PN 43975



www.airtechniques.com

Corporate Headquarters

1295 Walt Whitman Road | Melville, NY 11747-3062
Phone 800.247.8324 | Fax 888.247.8481

Western Facility

291 Bonnie Lane, Suite 101 | Corona, CA 92880-2804
Phone 800.247.8324 | Fax 951.898.7646

AIR
TECHNIQUES equipped for life™