

ScanX Swift View

by Air Techniques

The ScanX Swift View – the newest addition to the Air Techniques ScanX product portfolio – is a compact, single-slot phosphor storage plate (PSP) scanner digital radiography system. An ideal in-operatory solution or an adjunct to an existing direct digital sensor, the ScanX Swift View provides high-diagnostic image resolution (up to 20 lp/mm) with intuitive, time-saving digitization for all intraoral sizes (0-4). ScanX PSPs are as thin and flexible as traditional analog film, permitting easy placement for significantly increased patient comfort. And, greater patient comfort ensures fewer retakes, increasing productivity and profitability for the dental practice.



Dentists who are considering adding a ScanX Swift View should ask themselves several questions:

- Is my practice currently film-based, or digital?
- If I've already incorporated a direct-digital hard sensor, do I still find it difficult to capture images?"
- What percent of my patients are pediatric? What percent are geriatric?
- Do many of my patients have special needs or issues with gagging?
- Do I ever have the need to take a size 4 occlusal?

Some dentists may question the need to purchase a new direct-digital sensor. However, ScanX Swift View PSP technology should not be undervalued. Often, the best digital solution for a dental practice is a hybrid approach that includes both direct sensors and ScanX PSP technology. Dentists can capture 100 percent of their images with a hybrid digital solution, completely eliminating film and chemicals from the practice – a great solution for patients, the dental team and the practice.

Air Techniques offers a diverse product portfolio for today's dynamic market place, and the ScanX Swift View provides dentists with an innovative solution designed to address a variety of needs that may arise in their practice. ■

Editor's note: Sponsored by Air Techniques.

ONE SIZE DOES NOT FIT ALL!



ScanX PSPs are 30x thinner than rigid sensors.

*Which one is a better fit
for your patient?*

SCAN▶X®



What's all the chatter about? Visit
<https://www.airtechniques.com/scanx-chatter/>

For more information, to schedule a demo, or to place an order, contact your local dealer representative. To learn more, visit www.airtechniques.com

COMBINE SECTION 179 TAX INCENTIVES WITH PROMOTIONAL REBATES!
Help your business grow and deduct the FULL COST of Dental Equipment, Up to \$1,000,000! Please consult with your tax advisor for further information.



AIR
TECHNIQUES equipped for life®