



NUMBER OF EVALUATORS
13
COMBINED YEARS IN
PRACTICE
228

See what your colleagues have to say about hundreds of other evaluated products

View More Evaluations Online

www.dps.li/p/product-reviews

AIR TECHNIQUES

Monarch CleanStream

EVACUATION SYSTEM CLEANER

Monarch CleanStream Evacuation System Cleaner helps remove line buildup, increasing vacuum performance and suction.

A *Dental Product Shopper* evaluator with 31 years of experience that practices in Texas said of Monarch CleanStream, “It was extremely easy to use and it worked.” This evaluator rated all features as well as his overall satisfaction as excellent.

Monarch CleanStream is designed for the daily cleaning of dental evacuation systems, removing line buildup while increasing vacuum performance and suction. The unique, non-foaming formula cleans and deodorizes lines for all wet and dry vacuum systems, featuring an easy-to use dispenser system that makes mixing simple. Monarch CleanStream’s highly concentrated formula means a 2.5-L bottle yields 125 cleanings. It is also environmentally friendly and biodegradable.

This evaluation was conducted by a combination of current Monarch CleanStream users (3) and *Dental Product Shopper* evaluators (11).

Usage

Manufacturer Description: The Monarch CleanStream dispenser “automatically aspirates the solution into a fine mist for complete inner tubing wall coverage” without adding stress to the vacuum pump. The dispenser automatically dispenses the right amount of solution for each operatory, with no staff monitoring necessary. T



Evaluator Feedback: When asked to report their favorite feature of Monarch CleanStream, an evaluator with more than 20 years of experience in dentistry said, “I like that it was designed to be as adaptable to every operatory and customizable so as to make cleaning as easy as possible.”

An evaluator who rated his overall satisfaction with Monarch CleanStream as excellent reported, “[It is] easy to make the solution and the bottles are [easy] to attach to suction and let it do its thing.” Another evaluator said, “It was great to have ability to clean 2 operatories at once.”

One of the evaluators commented, “The large concentrate bottle could be improved. It gets messy,” and suggested that a “spout design” might be a possible improvement. However, this same evaluator still rated his overall satisfaction with Monarch CleanStream as excellent.

A few other evaluators mentioned that they liked the



It really works and it smells great.

Bruce Weitz, DDS
Lorain, OH

9 PREV
EVALS 35 YEARS
EXP



The container is awkward to use.

Annette Skowronski, DDS, FAGD
Shelby Township, MI

4 PREV
EVALS 28 YEARS
EXP




It was extremely easy to use and it worked.

Hal Board, DDS
Mesquite, TX


7 PREV
EVALS 31 YEARS
EXP

“shape of the dispenser” and said the product was “easy to mix.”



Vacuum Performance and Features

Manufacturer Description: CleanStream has excellent anti-foaming properties—the produced foam collapses immediately, increasing the longevity of pumps. Foam that passes through the separation tank can get into the suction pump, dry up, and cause the pump to fail. If foam gets into the plumbing, it can obstruct the flow of air, which creates the suction.



Evaluator Feedback: “It really works and it smells great,” said an evaluator who practices in Ohio, describing the performance of Monarch CleanStream. He added, “The staff likes it a lot.” Another evaluator said, “[The product] does the job it is intended to do.”

An evaluator with 9 years of experience said, “The anti-foaming properties are great,” and went on to rate all features and his overall satisfaction as excel-

lent. Another evaluator simply said, “It works,” and called Monarch CleanStream somewhat better than competing evacuation system cleaning products.

A dentist that practices in Wisconsin commented, “The product not leaving any colored residues makes it very attractive.”

Overall Satisfaction

Nine of the 13 evaluators rated their overall satisfaction with Monarch CleanStream as excellent or very good. One evaluator, who rated their overall satisfaction as excellent said he, “can’t

think of anything.” This evaluator said he would definitely recommend Monarch CleanStream to his colleagues and would probably purchase the product in the future.


FOR FREE INFORMATION:

800.420.8990 ext. XXXXX
www.dps.li/a/XXX-XXX
Circle XXX on the FREE info card following page XXX.



For more information scan the QR Code.

EVALUATION SNAPSHOT



FINAL SCORE
AVERAGE OF
SECTIONS A AND B

4.1

AIR TECHNIQUES, INC
Monarch CleanStream